



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 1/22/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 01/21/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:11
Out Date: 01/22/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB112870

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
M1112-07	S21-02-D	1	1.15	9.94	7.18	68.6	
M1112-08	S21-03-D	2	1.16	9.54	7.56	76.4	
M1112-09	S21-04-D	3	1.16	9.66	8.15	82.2	
M1113-20	S11-03-D	4	1.19	9.57	6.89	68.0	
M1113-21	S11-04-D	5	1.16	9.93	7.47	71.9	
M1113-22	S11-05-D	6	1.14	9.97	8.19	79.8	
M1113-23	S11-06-D	7	1.12	9.87	7.21	69.6	
M1113-24	S11-07-D	8	1.16	9.54	7.44	74.9	
M1113-25	S11-08-D	9	1.11	9.78	7.52	73.9	
M1113-26	S11-09-D	10	1.15	9.5	5.23	48.9	
M1113-27	DUP-53	11	1.15	9.96	8.41	82.4	
M1114-13	S22-03-D	12	1.16	9.93	8.25	80.8	
M1114-14	S22-04-D	13	1.15	9.76	7.71	76.2	
M1114-15	S22-05-D	14	1.19	9.77	7.42	72.6	
M1114-16	S22-06-D	15	1.1	9.62	6.53	63.7	
M1114-17	S22-07-D	16	1.19	9.82	8.6	85.9	
M1114-19	DUP-54	17	1.15	9.97	8.69	85.5	
M1115-09	S23-01-D	18	1.16	9.93	7.52	72.5	
M1115-10	S23-02-D	19	1.18	9.7	7.47	73.8	
M1115-11	S23-03-D	20	1.18	9.74	7.59	74.9	
M1115-12	S23-04-D	21	1.14	9.66	7.92	79.6	
M1185-01	KMA788P-END-1	40	1.00	2.00	2.00	100.0	PILC reel
M1185-02	KMA788P-END-2	41	1.00	2.00	2.00	100.0	PILC reel
M1186-01	KMB004E-END-1	42	1.00	2.00	2.00	100.0	PILC reel
M1186-02	KMB004E-END-2	43	1.00	2.00	2.00	100.0	PILC reel
M1187-01	K77029M-END-1	44	1.00	2.00	2.00	100.0	PILC reel
M1187-02	K77029M-END-2	45	1.00	2.00	2.00	100.0	PILC reel
M1187-03	ZZA182W-END-1	46	1.00	2.00	2.00	100.0	PILC reel



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 1/22/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 01/21/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:11
Out Date: 01/22/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB112870

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M1187-04	ZZA182W-END-2	47	1.00	2.00	2.00	100.0	PILC reel
M1187-05	216222-END-1	48	1.00	2.00	2.00	100.0	PILC reel
M1187-06	216222-END-2	49	1.00	2.00	2.00	100.0	PILC reel
M1187-07	UNKNOWN-MECHANIC-END-1	50	1.00	2.00	2.00	100.0	PILC reel
M1187-08	UNKNOWN-MECHANIC-END-2	51	1.00	2.00	2.00	100.0	PILC reel
M1188-01	TP-230	35	1.15	9.97	8.47	83.0	
M1188-02	TP-230	36	1.19	9.7	8.48	85.7	
M1188-04	TP-SWITCH-GEAR	37	1.13	9.58	8.23	84.0	
M1188-05	TP-SWITCH-GEAR	38	1.18	9.75	8.52	85.6	
M1188-05D	TP-SWITCH-GEARDUP	39	1.17	9.77	8.54	85.7	
M1191-01	0040-042-COMP-B	22	1.15	9.97	4.47	37.6	sediment sample
M1191-02	0041-043-COMP-C	23	1.1	9.68	4.4	38.5	sediment sample
M1191-03	0042-044-COMP-D	24	1.16	9.97	5.19	45.7	sediment sample
M1191-04	0043-045-COMP-E	25	1.17	9.64	5.22	47.8	sediment sample
M1191-05	0047-052-COMP-I	26	1.18	9.76	6.32	59.9	sediment sample
M1191-05D	0047-052-COMP-IDUP	27	1.19	9.78	6.34	60.0	sediment sample
M1193-01	OR-02-012121	33	1.18	9.68	8.8	89.6	
M1193-02	OR-02-012121-PEH-2	34	1.16	9.53	8.81	91.4	
M1197-01	VCP-EXC-SOIL	28	1.15	9.92	6.98	66.5	
M1198-01	TT-036-IDWSO-20210120	29	1.16	9.81	6.98	67.3	
M1198-02	TT-037-IDWSO-20210120	30	1.16	9.98	8.29	80.8	
M1200-01	012021TB01	31	1.00	2.00	2.00	100.0	CAULKING SAMPLE
M1200-02	012021TB02	32	1.00	2.00	2.00	100.0	CAULKING SAMPLE
M1206-01	FK-SF-03-009-A	52	1.12	9.73	8.77	88.9	
M1206-02	FK-SF-03-009-A	53	1.16	9.86	8.85	88.4	
M1206-03	FK-SF-03-009-B	54	1.13	9.65	8.68	88.6	
M1206-04	FK-SF-03-009-B	55	1.14	9.94	9.06	90.0	
M1207-01	FK-SF-03-010-A	56	1.17	9.6	8.58	87.9	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 1/22/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 01/21/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:11
Out Date: 01/22/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB112870

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M1207-02	FK-SF-03-010-A	57	1.1	9.7	8.77	89.2	
M1207-03	FK-SF-03-010-B	58	1.19	9.69	8.67	88.0	
M1207-04	FK-SF-03-010-B	59	1.14	9.81	8.85	88.9	
M1208-01	FK-SF-03-011-A	60	1.17	9.71	8.7	88.2	
M1208-02	FK-SF-03-011-A	61	1.18	9.98	9.14	90.5	
M1208-03	FK-SF-03-011-B	62	1.13	9.56	8.65	89.2	
M1208-04	FK-SF-03-011-B	63	1.19	9.84	8.86	88.7	
M1209-01	FK-PS-01-069	64	1.15	9.94	8.42	82.7	
M1209-01D	FK-PS-01-069DUP	65	1.15	9.96	8.45	82.9	
M1209-02	FK-PS-01-069	66	1.17	9.9	8.54	84.4	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$