



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 1/22/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 14:42
In Date: 01/21/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 104
Time OUT: 08:12
Out Date: 01/22/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB112878

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M1204-13	X2X96	1	1.00	2.00	2.00	100.0	P.E. SAMPLE
M1204-14	X2X96D	2	1.00	2.00	2.00	100.0	P.E. SAMPLE
M1204-15	X2X96S	3	1.00	2.00	2.00	100.0	P.E. SAMPLE
M1204-16	X2X97	4	1.00	2.00	2.00	100.0	P.E. SAMPLE
M1204-17	X2X97D	5	1.00	2.00	2.00	100.0	P.E. SAMPLE
M1204-18	X2X97S	6	1.00	2.00	2.00	100.0	P.E. SAMPLE
M1204-19	X2X99	7	1.00	2.00	2.00	100.0	P.E. SAMPLE
M1204-20	X2X99D	8	1.00	2.00	2.00	100.0	P.E. SAMPLE
M1204-21	X2X99S	9	1.00	2.00	2.00	100.0	P.E. SAMPLE
M1204-22	X2Y00	10	1.00	2.00	2.00	100.0	P.E. SAMPLE
M1204-23	X2Y00D	11	1.00	2.00	2.00	100.0	P.E. SAMPLE
M1204-24	X2Y00S	12	1.00	2.00	2.00	100.0	P.E. SAMPLE

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$