



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 1/25/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 17:15
In Date: 01/22/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:13
Out Date: 01/23/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB112890

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M1189-01	41108	1	1.00	2.00	2.00	100.0	OIL SAMPLE
M1189-02	41047-41048	2	1.18	9.7	9.08	92.7	
M1189-03	41041	3	1.00	2.00	2.00	100.0	OIL SAMPLE
M1189-04	41047	4	1.16	9.84	9.18	92.4	
M1189-05	41048	5	1.18	9.71	8.73	88.5	
M1213-01	IDW-BR2-SO-012121	6	1.15	9.95	7.88	76.5	
M1215-01	HR-02-012221	17	1.15	9.64	8.74	89.4	
M1215-02	HR-02-012221-EPH-2	18	1.12	9.61	8.74	89.8	
M1215-02D	HR-02-012221-EPH-2DUP	19	1.14	9.63	8.76	89.8	
M1216-01	PS-105	37	1.14	9.94	9.56	95.7	
M1216-02	PS-104	38	1.00	2.00	2.00	100.0	OILY DEBRIS
M1216-03	PS-102	39	1.00	2.00	2.00	100.0	OILY DEBRIS
M1217-01	1617	40	1.00	2.00	2.00	100.0	OILY DEBRIS
M1218-01	240	41	1.00	2.00	2.00	100.0	OILY DEBRIS
M1218-03	241	42	1.16	9.64	9.04	92.9	
M1220-01	SP1-A	43	1.14	9.71	8.7	88.2	
M1220-02	SP1-A-EPH-2	44	1.17	9.6	8.75	89.9	
M1220-03	SP1-B	45	1.16	9.7	9.08	92.7	
M1220-04	SP1-B-EPH-2	46	1.19	9.71	9.05	92.3	
M1220-05	SP1-C	47	1.18	9.74	8.69	87.7	
M1220-06	SP1-C-EPH-2	48	1.17	9.83	8.86	88.8	
M1221-05	SVOC-GPC-BLANK	7	1.00	2.00	2.00	100.0	
M1221-06	PEST-GPC-BLANK	8	1.00	2.00	2.00	100.0	
M1221-07	PEST-GPC-BLANK-SPIKE	9	1.00	2.00	2.00	100.0	
M1221-08	PCB-GPC-BLANK	10	1.00	2.00	2.00	100.0	
M1221-09	PCB-GPC-BLANK-SPIKE	11	1.00	2.00	2.00	100.0	
M1221-10	SVOC-GPC2-BLANK	12	1.00	2.00	2.00	100.0	
M1221-11	PEST-GPC2-BLANK	13	1.00	2.00	2.00	100.0	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 1/25/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 17:15
In Date: 01/22/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:13
Out Date: 01/23/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB112890

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
M1221-12	PEST-GPC2-BLANK-SPIKE	14	1.00	2.00	2.00	100.0	
M1221-13	PCB-GPC2-BLANK	15	1.00	2.00	2.00	100.0	
M1221-14	PCB-GPC2-BLANK-SPIKE	16	1.00	2.00	2.00	100.0	
M1224-01	BL-1	20	1.16	9.91	8.94	88.9	
M1224-02	BL-2	21	1.19	9.78	8.77	88.2	
M1224-03	BL-3	22	1.15	9.96	8.97	88.8	
M1224-04	BL-4	23	1.1	9.88	8.81	87.8	
M1224-05	BL-5	24	1.19	9.86	9.08	91.0	
M1224-06	BL-6	25	1.19	9.77	8.85	89.3	
M1224-07	BL-7	26	1.18	9.69	8.87	90.4	
M1224-08	BL-8	27	1.1	9.54	8.62	89.1	
M1225-01	S-1	28	1.16	9.54	8.64	89.3	
M1225-02	S-2	29	1.19	9.62	8.63	88.3	
M1225-03	S-3	30	1.17	9.74	8.87	89.8	
M1225-04	S-4	31	1.14	9.71	8.89	90.4	
M1225-05	S-5	32	1.15	9.54	8.85	91.8	
M1225-06	S-6	33	1.15	9.6	8.22	83.7	
M1225-07	S-7	34	1.15	9.54	8.74	90.5	
M1225-08	S-8	35	1.18	9.77	8.35	83.5	
M1225-08D	S-8DUP	36	1.15	9.9	8.37	82.5	
M1226-06	0046-051-COMP-H	49	1.16	9.98	7.96	77.1	
M1227-01	0065-3-2(0.0-0.5)	50	1.1	9.82	8.02	79.4	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$