



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 1/29/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 16:50
In Date: 01/28/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:02
Out Date: 01/29/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB112955

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M1267-01	1-5-TOTE-1-DRUM	1	1.00	2.00	2.00	100.0	OIL sample
M1267-02	1-3-DRUMS	2	1.00	2.00	2.00	100.0	oily rags
M1268-01	227	3	1.00	2.00	2.00	100.0	CONCRETE SAMPLE
M1269-01	CL-02-012821	11	1.16	9.94	9.25	92.1	
M1269-02	CL-02-012821-EPH-2	12	1.16	9.75	8.96	90.8	
M1270-01	RS-SP1	4	1.15	9.96	8.67	85.4	
M1270-02	RS-SP1-EPH-2	5	1.12	9.52	7.89	80.6	
M1270-03	RS-SP2	6	1.15	9.98	8.54	83.7	
M1270-04	RS-SP2-EPH-2	7	1.19	9.86	8.13	80.0	
M1270-05	RS-SP3	8	1.16	9.79	8.09	80.3	
M1270-06	RS-SP3-EPH-2	9	1.15	9.95	8.4	82.4	
M1270-06D	RS-SP3-EPH-2DUP	10	1.16	9.97	8.42	82.4	
M1271-01	RO-332	13	1.19	9.74	7.75	76.7	
M1271-02	RO-332	14	1.16	9.72	7.85	78.2	
M1271-03	RO-RBR-200035	15	1.1	9.78	8.55	85.8	
M1271-04	RO-RBR-200035	16	1.15	9.97	8.95	88.4	
M1279-01	S17-01-D	17	1.13	9.81	7.32	71.3	
M1279-02	S17-02-D	18	1.18	9.63	5.68	53.3	
M1279-03	S17-03-D	19	1.16	8.3	5.78	64.7	
M1279-04	S17-04-D	20	1.18	9.9	7.56	73.2	
M1279-05	S17-05-D	21	1.18	9.88	6.93	66.1	
M1279-06	S17-06-D	22	1.15	9.97	7.74	74.7	
M1279-07	S17-07-D	23	1.15	9.95	6.63	62.3	
M1280-01	S20-03-F	24	1.16	9.97	8.9	87.9	
M1280-02	S20-04-F	25	1.1	9.82	7.95	78.6	
M1280-03	S20-05-F	26	1.13	9.6	8.14	82.8	
M1280-04	DUP-60	27	1.19	9.58	7.3	72.8	
M1281-01	S21-02-E	28	1.13	9.9	7.46	72.2	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 1/29/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 16:50
In Date: 01/28/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:02
Out Date: 01/29/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB112955

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M1281-02	S21-03-E	29	1.16	9.97	8.18	79.7	
M1281-03	S21-04-E	30	1.16	9.98	7.65	73.6	
M1281-07	DUP-61	31	1.15	9.44	7.53	77.0	
M1282-01	S22-04-E	32	1.1	9.96	7.28	69.8	
M1282-02	S22-05-E	33	1.16	9.72	7.73	76.8	
M1282-03	S22-06-E	34	1.19	9.61	7.35	73.2	
M1282-04	DUP-62	35	1.12	9.68	7.13	70.2	
M1282-04D	DUP-62DUP	36	1.17	9.69	7.15	70.2	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$