



PERCENT SOLID

Supervisor: apatel
Analyst: Alexander
Date: 2/23/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 17:20
In Date: 02/22/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 105
Time OUT: 09:20
Out Date: 02/23/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB113216

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
M1460-01	HD-01-022221	1	1.15	9.74	8.06	80.4	
M1460-02	HD-01-022221-EPH-2	2	1.18	9.72	8.08	80.8	
M1461-01	31155	3	1.16	9.89	7.56	73.3	
M1462-01	CL-01-022221	4	1.16	9.88	8.28	81.7	
M1462-02	CL-01-022221-EPH-2	5	1.15	9.68	8.00	80.3	
M1464-01	2-22A	6	1.00	2.00	2.00	100.0	PLASTIC
M1464-02	2-22B	7	1.00	2.00	2.00	100.0	PLASTIC
M1464-03	2-22C	8	1.00	2.00	2.00	100.0	PLASTIC
M1464-03D	2-22CDUP	9	1.00	2.00	2.00	100.0	PLASTIC
M1464-04	2-22D	10	1.00	2.00	2.00	100.0	PLASTIC
M1464-05	2-22E	11	1.00	2.00	2.00	100.0	PLASTIC
M1465-01	KZZ094R-END-1	12	1.00	2.00	2.00	100.0	PILC REEL
M1465-02	KZZ094R-END-2	13	1.00	2.00	2.00	100.0	PILC REEL
M1465-03	HIB099Q-END-1	14	1.00	2.00	2.00	100.0	PILC REEL
M1465-04	HIB099Q-END-2	15	1.00	2.00	2.00	100.0	PILC REEL
M1466-01	JEC402U-END-1-1	16	1.00	2.00	2.00	100.0	PILC REEL
M1466-02	JEC402U-END-1-2	17	1.00	2.00	2.00	100.0	PILC REEL
M1466-03	JEC402U-END-2-1	18	1.00	2.00	2.00	100.0	PILC REEL
M1467-01	BC260098LOK-END-1	19	1.00	2.00	2.00	100.0	PILC REEL
M1467-02	BC260098LOK-END-2	20	1.00	2.00	2.00	100.0	PILC REEL
M1467-03	BC226716LOZ-END-1	21	1.00	2.00	2.00	100.0	PILC REEL
M1467-04	BC226716LOZ-END-2	22	1.00	2.00	2.00	100.0	PILC REEL

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$