



PERCENT SOLID

Supervisor: apatel
Analyst: Alexander
Date: 2/24/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 16:41
In Date: 02/23/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 106
Time OUT: 09:15
Out Date: 02/24/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB113233

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M1444-01	MGB2Y4	1	1.16	9.92	9.77	98.3	
M1444-02	MGB2Y5	2	1.17	9.91	9.62	96.7	
M1444-03	MGB2Y6	3	1.15	9.96	9.86	98.9	
M1444-04	MGB2Y7	4	1.18	9.6	9.44	98.1	
M1444-05	MGB2Y8	5	1.18	9.92	9.8	98.6	
M1444-06	MGB2Y9	6	1.17	9.85	9.77	99.1	
M1444-07	MGB2Y9D	7	1.17	9.85	9.77	99.1	
M1444-08	MGB2Y9S	8	1.17	9.85	9.77	99.1	
M1444-09	MGB2Z0	9	1.17	9.74	9.66	99.1	
M1444-10	MGB2Z1	10	1.18	9.63	9.51	98.6	
M1444-11	MGB2Z2	11	1.15	9.9	9.88	99.8	
M1444-12	MGB2Z3	12	1.00	2.00	2.00	100.0	P.E. SAMPLE

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$