



PERCENT SOLID

Supervisor: apatel
Analyst: Alexander
Date: 3/1/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 16:25
In Date: 02/24/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 106
Time OUT: 09:25
Out Date: 02/25/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB113249

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
M1475-01	WA-CON	1	1.15	9.75	9.05	91.9	
M1476-01	ML-CON	2	1.16	9.73	9.09	92.5	
M1477-01	72-11995	3	1.00	2.00	2.00	100.0	CONCRETE SAMPLE
M1477-02	72-11995-EPH-2	4	1.00	2.00	2.00	100.0	CONCRETE SAMPLE
M1478-01	SLABS-SLP	5	1.00	2.00	2.00	100.0	CONCRETE SAMPLE
M1480-01	JEH535U-END-1-1	6	1.00	2.00	2.00	100.0	PILC REEL
M1480-02	JEH535U-END-1-2	7	1.00	2.00	2.00	100.0	PILC REEL
M1480-03	KMH245E-END-1-1	8	1.00	2.00	2.00	100.0	PILC REEL
M1480-04	KMH245E-END-1-2	9	1.00	2.00	2.00	100.0	PILC REEL
M1480-05	KMH245E-END-2-1	10	1.00	2.00	2.00	100.0	PILC REEL
M1480-06	KMH245E-END-2-2	11	1.00	2.00	2.00	100.0	PILC REEL
M1480-07	GD1059I-END-1-1	12	1.00	2.00	2.00	100.0	PILC REEL
M1480-08	GD1059I-END-1-2	13	1.00	2.00	2.00	100.0	PILC REEL
M1480-09	HEC998P-END-1-1	14	1.00	2.00	2.00	100.0	PILC REEL
M1480-10	HEC998P-END-1-2	15	1.00	2.00	2.00	100.0	PILC REEL
M1480-11	HEC998P-END-2-1	16	1.00	2.00	2.00	100.0	PILC REEL
M1480-12	HEC998P-END-2-2	17	1.00	2.00	2.00	100.0	PILC REEL
M1480-13	KMH384W-END-1-1	18	1.00	2.00	2.00	100.0	PILC REEL
M1480-14	KMH384W-END-1-2	19	1.00	2.00	2.00	100.0	PILC REEL
M1480-15	KMH384W-END-2-1	20	1.00	2.00	2.00	100.0	PILC REEL
M1480-16	KMH384W-END-2-2	21	1.00	2.00	2.00	100.0	PILC REEL

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$