



PERCENT SOLID

Supervisor: apatel
Analyst: Alexander
Date: 3/1/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 11:30
In Date: 02/26/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:02
Out Date: 03/01/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB113279

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M1495-13	X3105	1	1.00	2.00	2.00	100.0	P.E. SAMPLE
M1495-14	X3105D	2	1.00	2.00	2.00	100.0	P.E. SAMPLE
M1495-15	X3105S	3	1.00	2.00	2.00	100.0	P.E. SAMPLE
M1495-16	X3106	4	1.00	2.00	2.00	100.0	P.E. SAMPLE
M1495-17	X3106D	5	1.00	2.00	2.00	100.0	P.E. SAMPLE
M1495-18	X3106S	6	1.00	2.00	2.00	100.0	P.E. SAMPLE
M1495-19	X3107	7	1.00	2.00	2.00	100.0	P.E. SAMPLE
M1495-20	X3107D	8	1.00	2.00	2.00	100.0	P.E. SAMPLE
M1495-21	X3107S	9	1.00	2.00	2.00	100.0	P.E. SAMPLE
M1495-22	X3108	10	1.00	2.00	2.00	100.0	P.E. SAMPLE
M1495-23	X3108D	11	1.00	2.00	2.00	100.0	P.E. SAMPLE
M1495-24	X3108S	12	1.00	2.00	2.00	100.0	P.E. SAMPLE

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$