



PERCENT SOLID

Supervisor: apatel
Analyst: Alexander
Date: 3/1/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 17:25
In Date: 02/26/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:00
Out Date: 03/01/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB113292

Table with 8 columns: Lab ID, Client SampleID, Dish #, Dish Wt (g) (A), Dish + Sample Wt (g) (B), Dish+Dry Sample Wt (g) (C), % Solid, Comments. Rows include samples like M1488-02, M1489-02, M1501-01, etc.



PERCENT SOLID

Supervisor: apatel
Analyst: Alexander
Date: 3/1/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 17:25
In Date: 02/26/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:00
Out Date: 03/01/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB113292

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$