



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 3/24/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 17:10
In Date: 03/23/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:22
Out Date: 03/24/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB113635

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M1730-01	TP-8A	1	1.14	9.93	8.17	80.0	
M1730-02	TP-8B	2	1.16	9.51	7.79	79.4	
M1769-01	SVE-03	3	1.00	2.00	2.00	100.0	OIL SAMPLE
M1771-01	COMP-SOIL-323	13	1.16	9.82	9.43	95.5	
M1771-02	COMP-SOIL-323-EPH-2	14	1.18	9.9	9.46	95.0	
M1772-01	DRUM-1-2-COMP	19	1.00	2.00	2.00	100.0	OILY DEBRIS
M1772-02	550026646	20	1.00	2.00	2.00	100.0	OIL SAMPLE
M1773-01	NB-08-032321	15	1.19	9.6	8.69	89.2	
M1773-02	NB-08-032321-EPH-2	16	1.19	9.58	8.75	90.1	
M1775-01	11979	26	1.15	9.98	3.02	21.2	GEL MATRIX
M1775-02	11979	27	1.14	9.85	3.37	25.6	gel matrix
M1776-01	110	28	1.00	2.00	2.00	100.0	OILY DEBRIS
M1776-03	106	29	1.00	2.00	2.00	100.0	OIL SAMPLE
M1777-01	244	30	1.15	9.97	7.61	73.2	
M1777-02	248-250	31	1.11	9.95	7.99	77.8	
M1779-01	SB-2B1- (3.0-3.5)	4	1.11	9.58	9.17	95.2	
M1779-02	SB-3B1- (3.0-3.5)	5	1.15	9.94	9.43	94.2	
M1779-04	SB-4B2- (4.0-4.5)	6	1.19	9.89	9.28	93.0	
M1779-05	SB-6B1- (3.0-3.5)	7	1.16	9.97	9.54	95.1	
M1779-06	SB-1B1- (3.0-3.5)	8	1.14	9.52	9.12	95.2	
M1779-07	SB-1B2- (3.5-4.0)	9	1.16	9.96	9.55	95.3	
M1779-08	DUP032221	10	1.18	9.95	9.56	95.6	
M1780-01	TT-41-RW6-IDWSO-20210322	11	1.18	9.58	9.12	94.5	
M1780-01D	TT-41-RW6-IDWSO-20210322DUP	12	1.19	9.6	9.14	94.5	
M1786-01	SP-1	17	1.16	9.7	8.91	90.7	
M1786-02	SP-1A	18	1.16	9.93	9.1	90.5	
M1789-01	PIPE-02	21	1.11	9.88	8.97	89.6	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 3/24/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 17:10
In Date: 03/23/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:22
Out Date: 03/24/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB113635

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M1789-02	M1789-01MS	22	1.11	9.88	8.97	89.6	
M1789-03	M1789-01MSD	23	1.11	9.88	8.97	89.6	
M1789-04	PIPE-01	24	1.13	9.81	8.79	88.2	
M1789-05	PIPE-019	25	1.17	9.75	8.68	87.5	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$