



PERCENT SOLID

Supervisor: apatel  
Analyst: jignesh  
Date: 4/7/2021

OVENTEMP IN Celsius(°C): 107  
Time IN: 14:45  
In Date: 04/01/2021  
Weight Check 1.0g: 1.00  
Weight Check 10g: 10.00  
OvenID: M Oven-1

OVENTEMP OUT Celsius(°C): 103  
Time OUT: 08:00  
Out Date: 04/02/2021  
Weight Check 1.0g: 1.00  
Weight Check 10g: 10.00  
BalanceID: M SC-4  
Thermometer ID: %solids-oven

QC:LB113784

| Lab ID   | Client SampleID | Dish # | Dish Wt (g) (A) | Dish + Sample Wt (g) (B) | Dish+Dry Sample Wt (g) (C) | % Solid | Comments    |
|----------|-----------------|--------|-----------------|--------------------------|----------------------------|---------|-------------|
| M1904-13 | X3248           | 1      | 1.00            | 2.00                     | 2.00                       | 100.0   | P.E. SAMPLE |
| M1904-14 | X3248D          | 2      | 1.00            | 2.00                     | 2.00                       | 100.0   | P.E. SAMPLE |
| M1904-15 | X3248S          | 3      | 1.00            | 2.00                     | 2.00                       | 100.0   | P.E. SAMPLE |
| M1904-16 | X3249           | 4      | 1.00            | 2.00                     | 2.00                       | 100.0   | P.E. SAMPLE |
| M1904-17 | X3249D          | 5      | 1.00            | 2.00                     | 2.00                       | 100.0   | P.E. SAMPLE |
| M1904-18 | X3249S          | 6      | 1.00            | 2.00                     | 2.00                       | 100.0   | P.E. SAMPLE |
| M1904-19 | X3250           | 7      | 1.00            | 2.00                     | 2.00                       | 100.0   | P.E. SAMPLE |
| M1904-20 | X3250D          | 8      | 1.00            | 2.00                     | 2.00                       | 100.0   | P.E. SAMPLE |
| M1904-21 | X3250S          | 9      | 1.00            | 2.00                     | 2.00                       | 100.0   | P.E. SAMPLE |
| M1904-22 | X3251           | 10     | 1.00            | 2.00                     | 2.00                       | 100.0   | P.E. SAMPLE |
| M1904-23 | X3251D          | 11     | 1.00            | 2.00                     | 2.00                       | 100.0   | P.E. SAMPLE |
| M1904-24 | X3251S          | 12     | 1.00            | 2.00                     | 2.00                       | 100.0   | P.E. SAMPLE |

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$