



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 4/15/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 17:25
In Date: 04/14/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:20
Out Date: 04/15/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB113973

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M2028-01	HC-1-OIL	38	1.00	2.00	2.00	100.0	OILY SAMPLE
M2028-03	HC-2-7-OILY-DEBRIS	39	1.00	2.00	2.00	100.0	OILY DEBRIS
M2028-05	HC-8	40	1.00	2.00	2.00	100.0	oily-stone
M2028-07	HC-9	41	1.00	2.00	2.00	100.0	CABLE SAMPLE
M2029-01	OR-01-041421	30	1.12	9.74	8.83	89.4	
M2029-02	OR-01-041421-EPH-2	31	1.12	9.66	8.94	91.6	
M2030-01	04-14-21-A	28	1.00	2.00	2.00	100.0	OILY DEBRIS
M2030-03	ETGI-304	29	1.00	2.00	2.00	100.0	PLASTIC
M2034-02	1979	42	1.17	9.81	8.3	82.5	
M2034-03	1979	43	1.13	9.75	7.94	79.0	
M2036-01	RT-80-W-ROCKAWAY-A1	1	1.12	9.88	8.62	85.6	
M2036-02	RT-80-W-ROCKAWAY-A1	2	1.12	9.72	8.19	82.2	
M2036-03	RT-80-W-ROCKAWAY-A2	3	1.14	9.94	8.6	84.8	
M2036-04	RT-80-W-ROCKAWAY-A2	4	1.16	9.83	8.45	84.1	
M2036-05	RT-80-W-ROCKAWAY-B1	5	1.18	9.59	8.26	84.2	
M2036-05D	RT-80-W-ROCKAWAY-B1DUP	6	1.16	9.61	8.29	84.4	
M2036-06	RT-80-W-ROCKAWAY-B1	7	1.18	9.68	8.33	84.1	
M2036-07	RT-80-W-ROCKAWAY-B2	8	1.19	9.73	8.29	83.1	
M2036-08	RT-80-W-ROCKAWAY-B2	9	1.15	9.52	8.24	84.7	
M2040-01	FILL-A-3	32	1.13	6.91	5.83	81.3	
M2040-03	FILL-A	33	1.16	9.84	8.11	80.1	
M2040-05	FILL-B-3	34	1.12	5.98	5.51	90.3	
M2040-07	FILL-B	35	1.18	9.79	8.16	81.1	
M2040-09	FILL-C-2	36	1.18	7.15	6.11	82.6	
M2040-11	FILL-C	37	1.13	9.96	8.72	86.0	
M2041-01	S-2A	10	1.15	9.79	8.47	84.7	
M2041-02	S-2B	11	1.16	9.76	7.9	78.4	
M2041-03	S-2C	12	1.1	9.66	7.97	80.3	



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 4/15/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 17:25
In Date: 04/14/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:20
Out Date: 04/15/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB113973

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M2041-04	S-2D	13	1.15	9.95	8.78	86.7	
M2042-01	FS-SB12A (7-8)	22	1.13	9.77	9.26	94.1	
M2042-02	M2042-01MS	23	1.13	9.77	9.26	94.1	
M2042-03	M2042-01MSD	24	1.13	9.77	9.26	94.1	
M2042-04	FS-SB12A (41-42)	25	1.14	9.71	9.48	97.3	
M2042-05	FS-SB12A (55-56)	26	1.13	9.86	8.98	89.9	
M2042-06	FS-SB12A (61-62)	27	1.16	9.98	9.05	89.5	
M2044-01	TR-05-041421	44	1.12	9.71	8.76	88.9	
M2044-02	TR-05-041421-EPH-2	45	1.13	9.86	8.84	88.3	
M2044-02D	TR-05-041421-EPH-2DUP	46	1.19	9.88	8.86	88.3	
M2046-01	WC-4	14	1.15	9.92	9.33	93.3	
M2046-02	WC-5	15	1.16	9.81	8.94	89.9	
M2047-01	130-SB2W- (9.5-10)	16	1.12	9.53	8.57	88.6	
M2047-02	130-SB2N- (9.5-10)	17	1.16	9.64	8.86	90.8	
M2047-03	130-SB2S- (9.5-10)	18	1.12	9.94	9.34	93.2	
M2047-04	130-SB2E- (9.5-10)	19	1.15	9.52	8.74	90.7	
M2047-05	130-SB3- (1.5-2)	20	1.12	9.98	7.93	76.9	
M2047-06	124-SB-1	21	1.17	9.98	8.48	83.0	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$