



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 4/16/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 17:05
In Date: 04/15/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:15
Out Date: 04/16/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID:

QC:LB113988

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M2043-01	FS-WC-SOIL	1	1.17	9.58	8.5	87.2	
M2043-01D	FS-WC-SOILDUP	2	1.18	9.98	8.52	83.4	
M2048-01	SB-01-0-2.0	3	1.17	9.69	8.52	86.3	
M2048-02	SB-01-COMP	4	1.16	9.94	8.46	83.1	
M2048-03	SB-01-14.5-15.0	5	1.17	9.57	8.23	84.0	
M2048-04	SB-02-0-2.0	6	1.17	9.87	8.24	81.3	
M2048-05	SB-02-COMP	7	1.18	9.89	8.71	86.5	
M2048-06	SB-02-13.5-14.0	8	1.17	9.52	8.41	86.7	
M2048-07	SB-03-0-2.0	9	1.16	9.59	8.55	87.7	
M2048-08	SB-03-COMP	10	1.17	9.51	8.24	84.8	
M2048-09	SB-03-14.5-15.0	11	1.17	9.52	8.77	91.0	
M2048-10	SB-04-0-2.0	12	1.16	9.7	8.43	85.1	
M2048-11	SB-04-COMP	13	1.18	9.98	8.73	85.8	
M2048-12	SB-04-14.5-15.0	14	1.17	9.68	8.51	86.3	
M2048-13	SB-05-0-2.0	15	1.16	9.61	9.14	94.4	
M2048-14	SB-05-COMP	16	1.18	9.55	8.1	82.7	
M2048-15	SB-05-14.5-15.0	17	1.18	9.79	9.06	91.5	
M2048-16	SB-06-0-2.0	18	1.17	9.85	8.8	87.9	
M2048-17	SB-06-COMP	19	1.17	9.59	8.53	87.4	
M2048-18	SB-06-8.5-9.0	20	1.16	9.82	8.36	83.1	
M2048-19	SB-07-0-2.0	21	1.18	9.93	9.32	93.0	
M2048-20	SB-07-COMP	22	1.18	9.75	9.19	93.5	
M2048-20D	SB-07-COMPDUP	23	1.19	9.74	9.2	93.7	
M2049-01	SB-07-9.5-10.0	24	1.18	9.81	8.49	84.7	
M2049-02	SB-08-0-2.0	25	1.17	9.87	9.26	93.0	
M2049-03	SB-08-COMP	26	1.18	9.91	9.04	90.0	
M2049-04	SB-08-7.5-8.0	27	1.18	9.87	8.82	87.9	
M2049-05	SB-08-7.5-8.0-DUP	28	1.17	9.81	8.94	89.9	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 4/16/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 17:05
In Date: 04/15/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:15
Out Date: 04/16/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID:

QC:LB113988

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M2049-06	SB-09-0-2.0	29	1.19	9.92	9.24	92.2	
M2049-07	SB-09-COMP	30	1.17	9.86	7.58	73.8	
M2049-08	SB-09-5.5-6.0	31	1.18	9.85	9.24	93.0	
M2049-10	SB-10-0-2.0	32	1.18	9.79	9.13	92.3	
M2049-11	SB-10-COMP	33	1.19	9.9	9.04	90.1	
M2049-12	SB-10-8.5-9.0	34	1.17	9.52	8.54	88.3	
M2049-13	SB-11-0-2.0	35	1.17	9.55	9.06	94.2	
M2049-14	SB-11-COMP	36	1.18	9.51	8.81	91.6	
M2049-15	SB-11-8.5-9.0	37	1.16	9.93	9.42	94.2	
M2049-16	SB-12-0-2.0	38	1.17	9.9	9.14	91.3	
M2049-17	SB-12-COMP	39	1.17	9.76	8.86	89.5	
M2049-18	SB-12-13.5-14.0	40	1.17	9.99	9.44	93.8	
M2050-01	TP-7	41	1.17	9.57	7.81	79.0	
M2050-03	TP-8	42	1.17	9.67	7.98	80.1	
M2050-05	TP-9	43	1.16	9.68	8.17	82.3	
M2050-05D	TP-9DUP	44	1.16	9.52	8.19	84.1	
M2050-07	TP-5	45	1.18	9.74	8.76	88.6	
M2050-09	TP-4	46	1.15	9.55	7.98	81.3	
M2050-11	TP-3	47	1.17	9.88	8.04	78.9	
M2050-13	TP-2	48	1.17	9.65	7.73	77.4	
M2050-15	TP-1	49	1.17	9.69	8.46	85.6	
M2050-17	TP-6	50	1.16	9.66	8.26	83.5	
M2051-01	0124027	51	1.16	6.62	5.43	78.2	
M2051-02	0124027-EPH-2	52	1.16	9.83	7.36	71.5	
M2051-03	0124035	53	1.16	9.98	6.94	65.5	
M2051-04	0124035-EPH-2	54	1.16	9.71	6.96	67.8	
M2051-05	0124055	55	1.17	9.91	7.08	67.6	
M2051-06	0124055-EPH-2	56	1.17	9.76	7.23	70.5	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 4/16/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 17:05
In Date: 04/15/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:15
Out Date: 04/16/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID:

QC:LB113988

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M2051-07	0124019	57	1.18	4.78	3.75	71.4	
M2051-08	0124019-EPH-2	58	1.16	9.66	7.03	69.1	
M2052-01	OK-02-041521	59	1.17	9.57	9.18	95.4	
M2052-02	OK-02-041521-EPH-2	60	1.18	9.85	9.4	94.8	
M2053-01	PADS-2	61	1.00	2.00	2.00	100.0	PADS
M2054-01	LONGEST-PIPE-1-2	62	1.00	2.00	2.00	100.0	WIPE SAMPLE
M2054-02	LONGEST-PIPE-2-2	63	1.00	2.00	2.00	100.0	WIPE SAMPLE

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$