



PERCENT SOLID

Supervisor: ketankumar

Analyst: JIGNESH

Date: 4/17/2021

OVENTEMP IN Celsius(°C): 107

Time IN: 17:20

In Date: 04/16/2021

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103

Time OUT: 08:11

Out Date: 04/17/2021

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

BalanceID: M SC-4

Thermometer ID: %SOLIDS-OVEN

QC:LB114016

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
M2063-01	TP1-0-2	1	1.15	9.96	8.84	87.3	
M2063-02	TP1-0-2-EPH-2	2	1.16	9.65	8.55	87.0	
M2063-03	TP1-2-4	3	1.15	9.67	8.72	88.8	
M2063-04	TP1-2-4-EPH-2	4	1.14	9.95	9.02	89.4	
M2063-05	TP2-0-2	5	1.18	9.58	8.44	86.4	
M2063-06	TP2-0-2-EPH-2	6	1.16	9.87	8.7	86.6	
M2063-07	TP2-2-4	7	1.18	9.64	8.66	88.4	
M2063-08	TP2-2-4-EPH-2	8	1.12	9.89	8.7	86.4	
M2063-09	TP3-0-2	9	1.17	9.73	8.9	90.3	
M2063-10	TP3-0-2-EPH-2	10	1.14	9.9	8.86	88.1	
M2063-11	TP3-2-4	11	1.16	9.93	9.02	89.6	
M2063-12	TP3-2-4-EPH-2	12	1.13	9.87	8.97	89.7	
M2063-12D	TP3-2-4-EPH-2DUP	13	1.14	9.89	8.99	89.7	
M2064-15	SVOC-GCP-BLANK	14	1.00	2.00	2.00	100.0	
M2064-16	PEST-GPC-BLANK	15	1.00	2.00	2.00	100.0	
M2064-17	PEST-GPC-BLANK-SPIKE	16	1.00	2.00	2.00	100.0	
M2064-18	PCB-GPC-BLANK	17	1.00	2.00	2.00	100.0	
M2064-19	PCB-GPC-BLANK-SPIKE	18	1.00	2.00	2.00	100.0	
M2064-20	SVOC-GCP2-BLANK	19	1.00	2.00	2.00	100.0	
M2064-21	PEST-GCP2-BLANK	20	1.00	2.00	2.00	100.0	
M2064-22	PEST-GCP2-BLANK-SPIKE	21	1.00	2.00	2.00	100.0	
M2064-23	PCB-GCP2-BLANK	22	1.00	2.00	2.00	100.0	
M2064-24	PCB-GCP2-BLANK-SPIKE	23	1.00	2.00	2.00	100.0	
M2068-01	PCB-041421-01	68	1.00	2.00	2.00	100.0	CAULKING
M2068-02	PCB-041421-02	69	1.00	2.00	2.00	100.0	CAULKING
M2068-03	PCB-041421-03	70	1.00	2.00	2.00	100.0	CAULKING
M2068-04	PCB-041421-04	71	1.00	2.00	2.00	100.0	CAULKING
M2068-05	PCB-041421-05	72	1.00	2.00	2.00	100.0	CAULKING



PERCENT SOLID

Supervisor: ketankumar

Analyst: JIGNESH

Date: 4/17/2021

OVENTEMP IN Celsius(°C): 107

Time IN: 17:20

In Date: 04/16/2021

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103

Time OUT: 08:11

Out Date: 04/17/2021

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

BalanceID: M SC-4

Thermometer ID: %SOLIDS-OVEN

QC:LB114016

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M2070-01	JGD-3150-D3550-A	24	1.16	9.76	3.75	30.1	SLUDGE SAMPLE
M2070-03	JGD-3150-D3550-A	25	1.13	9.84	3.94	32.3	SLUDGE SAMPLE
M2070-05	JGD-3150-D3550-B	26	1.15	9.95	4.07	33.2	SLUDGE SAMPLE
M2070-07	JGD-3150-D3550-B	27	1.17	9.97	4.25	35.0	SLUDGE SAMPLE
M2074-01	VNPH2-B-1	28	1.18	9.74	7.37	72.3	
M2074-02	VNPH2-B-1D	29	1.16	9.98	6.88	64.9	
M2074-03	VNPH2-B-2	30	1.14	9.94	8.7	85.9	
M2074-04	VNPH2-B-2D	31	1.19	9.56	6.66	65.4	
M2074-05	VNPH2-B-3	32	1.14	9.94	9.1	90.5	
M2074-06	VNPH2-B-3D	33	1.11	9.88	6.23	58.4	
M2074-07	VNPH2-B-4	34	1.15	9.96	8.85	87.4	
M2074-08	VNPH2-B-5	35	1.13	9.51	8.15	83.8	
M2074-09	VNPH2-B-6	36	1.13	9.74	8.09	80.8	
M2074-10	VNPH2-B-6D	37	1.14	9.97	8.28	80.9	
M2074-11	VNPH2-B-7	38	1.14	9.95	8.44	82.9	
M2074-12	VNPH2-B-8	39	1.16	9.85	7.87	77.2	
M2074-13	VNPH2-B-9	40	1.15	9.52	8.11	83.2	
M2074-14	VNPH2-B-10	41	1.16	9.92	8.73	86.4	
M2074-14D	VNPH2-B-10DUP	42	1.14	9.94	8.75	86.5	
M2074-15	VNPH2-B-10D	43	1.19	9.64	6.4	61.7	
M2074-16	VNPH2-B-11	44	1.16	9.93	8.83	87.5	
M2074-17	VNPH2-B-12	45	1.11	9.95	8.43	82.8	
M2075-01	TP-1-S	46	1.14	9.95	8.88	87.9	
M2075-03	TP-1-S1	47	1.14	9.9	9.02	90.0	
M2075-04	TP-1-D	48	1.19	9.98	9.09	89.9	
M2075-06	TP-1-D1	49	1.16	9.78	8.99	90.8	
M2075-07	TP-2-S	50	1.16	9.77	8.74	88.0	
M2075-09	TP-2-S1	51	1.15	9.95	8.79	86.8	



PERCENT SOLID

Supervisor: ketankumar

Analyst: JIGNESH

Date: 4/17/2021

OVENTEMP IN Celsius(°C): 107

Time IN: 17:20

In Date: 04/16/2021

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103

Time OUT: 08:11

Out Date: 04/17/2021

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

BalanceID: M SC-4

Thermometer ID: %SOLIDS-OVEN

QC:LB114016

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M2075-10	TP-2-D	52	1.12	9.84	9.06	91.1	
M2075-12	TP-2-D1	53	1.11	9.57	8.9	92.1	
M2075-13	TP-4-S	54	1.13	9.74	8.38	84.2	
M2075-15	TP-4-S1	55	1.14	9.94	8.97	89.0	
M2075-16	TP-4-D	56	1.14	9.95	9.18	91.3	
M2075-18	TP-4-D1	57	1.13	9.92	9.00	89.5	
M2075-19	TP-6-S	58	1.14	9.97	8.54	83.8	
M2075-21	TP-6-S1	59	1.17	9.68	8.2	82.6	
M2075-22	TP-6-D	60	1.14	9.9	9.02	90.0	
M2075-24	TP-6-D1	61	1.15	9.97	9.34	92.9	
M2075-25	TP-8-S	62	1.13	9.75	9.00	91.3	
M2075-27	TP-8S1	63	1.16	9.96	9.22	91.6	
M2075-28	TP-8-D	64	1.16	9.54	8.83	91.5	
M2075-30	TP-8-D1	65	1.13	9.7	8.95	91.2	
M2075-30D	TP-8-D1DUP	66	1.16	9.72	8.97	91.2	
M2076-01	SB-13-0-2.0	73	1.16	9.87	9.34	93.9	
M2076-02	SB-13-COMP	74	1.16	9.82	9.08	91.5	
M2076-03	SB-13-14.5-15.0	75	1.13	9.96	9.05	89.7	
M2076-04	SB-14-0-2.0	76	1.16	9.55	9.19	95.7	
M2076-05	SB-14-COMP	77	1.15	9.94	9.24	92.0	
M2076-06	SB-14-14.5-15.0	78	1.13	9.76	9.04	91.7	
M2076-07	SB-15-0-2.0	79	1.16	9.5	8.42	87.1	
M2076-08	SB-15-COMP	80	1.19	9.69	8.82	89.8	
M2076-09	SB-15-9.5-10.0	81	1.17	9.61	8.83	90.8	
M2076-10	SB-16-0-2.0	82	1.2	9.57	8.61	88.5	
M2076-11	SB-16-COMP	83	1.19	9.77	8.86	89.4	
M2076-12	SB-16-9.5-10.0	84	1.13	9.58	8.95	92.5	
M2076-13	SB-17-0-2.0	85	1.1	9.69	8.25	83.2	



PERCENT SOLID

Supervisor: ketankumar

Analyst: JIGNESH

Date: 4/17/2021

OVENTEMP IN Celsius(°C): 107

Time IN: 17:20

In Date: 04/16/2021

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103

Time OUT: 08:11

Out Date: 04/17/2021

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

BalanceID: M SC-4

Thermometer ID: %SOLIDS-OVEN

QC:LB114016

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M2076-14	SB-17-COMP	86	1.18	9.74	8.76	88.6	
M2076-15	SB-17-14.5-15.0	87	1.19	9.84	9.06	91.0	
M2076-16	SB-18-0-2.0	88	1.18	9.54	8.69	89.8	
M2076-17	SB-18-COMP	89	1.15	9.95	9.24	91.9	
M2076-18	SB-18-14.5-15.0	90	1.14	9.97	9.2	91.3	
M2076-19	SB-19-0-2.0	91	1.16	9.81	8.84	88.8	
M2076-20	SB-19-COMP	92	1.16	9.61	8.55	87.5	
M2077-01	SB-19-14.5-15.0	93	1.16	9.97	8.62	84.7	
M2077-02	SB-20-0-2.0	94	1.18	9.98	9.00	88.9	
M2077-03	SB-20-COMP	95	1.13	9.97	8.67	85.3	
M2077-04	SB-20-14.5-15.0	96	1.18	9.63	8.55	87.2	
M2077-05	SB-20-14.5-15.0-DUP	97	1.13	9.85	8.47	84.2	
M2077-06	SB-21-0-2.0	98	1.15	9.96	9.34	93.0	
M2077-07	SB-21-COMP	99	1.14	9.9	8.69	86.2	
M2077-08	SB-21-14.5-15.0	100	1.11	9.74	8.8	89.1	
M2077-08D	SB-21-14.5-15.0DUP	101	1.17	9.76	8.82	89.1	
M2080-01	2516	67	1.00	2.00	2.00	100.0	CONCRETE SAMPLE

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$