



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 4/21/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 04/20/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:19
Out Date: 04/21/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB114071

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M2098-01	HR-01-042021	1	1.16	9.55	8.51	87.6	
M2098-02	HR-01-042021-EPH-2	2	1.16	9.84	8.89	89.1	
M2098-03	HR-04-042021	3	1.12	9.82	8.74	87.6	
M2098-04	HR-04-042021-EPH-2	4	1.13	9.79	9.42	95.7	
M2099-01	TP-1	73	1.14	9.97	9.04	89.5	
M2099-02	TP-1-EPH-2	74	1.17	9.74	8.58	86.5	
M2099-03	TP-2	75	1.15	9.98	9.09	89.9	
M2099-04	TP-2-EPH-2	76	1.18	9.79	8.97	90.5	
M2099-05	TP-3	77	1.18	9.66	8.51	86.4	
M2099-06	TP-3-EPH-2	78	1.11	9.77	8.12	80.9	
M2099-07	TP-4	79	1.15	9.97	9.3	92.4	
M2099-08	TP-4-EPH-2	80	1.16	9.85	8.74	87.2	
M2099-09	TP-5	81	1.14	9.94	9.16	91.1	
M2099-10	TP-5-EPH-2	82	1.14	9.51	8.49	87.8	
M2100-01	361	5	1.14	9.96	8.6	84.6	
M2100-02	361	6	1.16	9.77	8.73	87.9	
M2101-01	31042-43	7	1.19	9.77	8.03	79.7	
M2101-02	31045	8	1.16	9.79	8.15	81.0	
M2101-03	30837	9	1.1	9.78	8.82	88.9	
M2102-01	RO-124037	83	1.17	9.72	7.74	76.8	
M2102-03	RB37708RT	84	1.19	9.72	7.49	73.9	
M2102-05	RB44882RT	85	1.15	9.98	8.57	84.0	
M2103-01	RO-124030	86	1.15	9.55	7.73	78.3	
M2103-01D	RO-124030DUP	87	1.18	9.58	7.75	78.2	
M2106-01	SB-22-0-2.0	10	1.12	9.87	8.89	88.8	
M2106-02	SB-22-COMP	11	1.16	9.66	8.74	89.2	
M2106-03	SB-22-14.5-15.0	12	1.19	9.81	8.67	86.8	
M2106-04	SB-23-0-2.0	13	1.1	9.81	8.9	89.6	



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 4/21/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 04/20/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:19
Out Date: 04/21/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB114071

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M2106-05	SB-23-COMP	14	1.12	9.71	8.65	87.7	
M2106-06	SB-23-14.5-15.0	15	1.13	9.63	8.78	90.0	
M2106-07	SB-24-0-2.0	16	1.15	9.98	9.00	88.9	
M2106-08	SB-24-COMP	17	1.19	9.66	8.32	84.2	
M2106-09	SB-24-14.5-15.0	18	1.15	9.97	8.46	82.9	
M2106-10	SB-25-0-2.0	19	1.16	9.69	8.59	87.1	
M2106-11	SB-25-COMP	20	1.15	9.94	8.8	87.0	
M2106-11D	SB-25-COMPDUP	21	1.16	9.96	8.82	87.0	
M2106-12	SB-25-9.5-10.0	22	1.13	9.71	8.00	80.1	
M2106-13	SB-26-0-2.0	23	1.16	9.88	8.94	89.2	
M2106-14	SB-26-COMP	24	1.12	9.54	8.48	87.4	
M2106-15	SB-26-7.5-8.0	25	1.11	9.82	8.32	82.8	
M2106-16	SB-27-0-2.0	26	1.14	9.85	8.5	84.5	
M2106-17	SB-27-COMP	27	1.15	9.97	9.45	94.1	
M2106-18	SB-27-9.5-10.0	28	1.13	9.85	8.19	81.0	
M2106-19	SB-28-0-2.0	29	1.1	9.88	9.06	90.7	
M2106-20	SB-28-COMP	30	1.15	9.6	8.68	89.1	
M2106-21	SB-28-9.5-10.0	31	1.16	9.8	8.46	84.5	
M2107-01	SB-29-0-2.0	32	1.14	9.71	8.91	90.7	
M2107-02	SB-29-COMP	33	1.16	9.55	8.57	88.3	
M2107-03	SB-29-9.5-10.0	34	1.19	9.66	8.73	89.0	
M2107-04	SB-30-0-2.0	35	1.18	9.98	9.27	91.9	
M2107-05	SB-30-COMP	36	1.14	9.93	8.97	89.1	
M2107-06	SB-30-9.5-10.0	37	1.1	9.89	9.19	92.0	
M2107-07	SB-30-9.5-10.0-DUP	38	1.15	9.88	8.88	88.5	
M2107-08	SB-31-0-2.0	39	1.17	9.6	8.55	87.5	
M2107-09	SB-31-COMP	40	1.17	9.77	8.62	86.6	
M2107-09D	SB-31-COMPDUP	41	1.15	9.79	8.64	86.7	



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 4/21/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 04/20/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:19
Out Date: 04/21/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB114071

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M2107-10	SB-31-10.5-11.0	42	1.15	9.98	9.09	89.9	
M2107-11	SB-32-0-2.0	43	1.13	9.72	8.21	82.4	
M2107-12	SB-32-COMP	44	1.14	9.94	8.74	86.4	
M2107-13	SB-32-10.5-11.0	45	1.12	9.58	8.45	86.6	
M2107-14	SB-33-0-2.0	46	1.13	9.82	8.44	84.1	
M2107-15	SB-33-COMP	47	1.14	9.94	8.06	78.6	
M2107-16	SB-33-10.5-11.0	48	1.18	9.98	8.58	84.1	
M2107-17	SB-34-0-2.0	49	1.16	9.74	8.6	86.7	
M2107-18	SB-34-COMP	50	1.15	9.95	8.53	83.9	
M2107-19	SB-34-9.5-10.0	51	1.18	9.92	8.02	78.3	
M2107-20	SB-35-0-2.0	52	1.11	9.65	9.06	93.1	
M2107-21	SB-35-COMP	53	1.18	9.67	8.99	92.0	
M2107-22	SB-35-10.5-11.0	54	1.13	9.64	7.4	73.7	
M2109-01	PE-1	55	1.18	9.7	8.3	83.6	
M2109-02	PE-2	56	1.14	9.98	8.83	87.0	
M2109-03	PE-3	57	1.13	9.58	8.1	82.5	
M2109-04	PE-4	58	1.12	9.9	8.4	82.9	
M2109-05	PE-5	59	1.17	9.9	8.63	85.5	
M2110-01	S-5	60	1.17	9.59	8.94	92.3	
M2110-02	S-6	61	1.18	9.7	9.00	91.8	
M2110-03	S-7	62	1.18	9.64	8.96	92.0	
M2110-04	S-8	63	1.1	9.9	9.21	92.2	
M2111-01	S-1	64	1.14	9.97	9.42	93.8	
M2111-02	S-2	65	1.12	9.91	9.31	93.2	
M2111-03	S-3	66	1.19	9.73	9.02	91.7	
M2111-04	S-4	67	1.15	9.55	8.8	91.1	
M2111-05	S-9	68	1.15	9.98	8.88	87.5	
M2111-06	S-10	69	1.13	9.91	9.33	93.4	



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 4/21/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 04/20/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:19
Out Date: 04/21/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB114071

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M2111-07	S-11	70	1.18	9.6	9.11	94.2	
M2112-01	M-1	71	1.15	9.96	8.85	87.4	
M2112-02	M-2	72	1.19	9.77	9.36	95.2	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$