



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 4/23/2021

OVENTEMP IN Celsius(°C): 106
Time IN: 17:00
In Date: 04/22/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:22
Out Date: 04/23/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB114113

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M2132-01	TP-1-COMP	1	1.15	9.85	8.99	90.1	
M2132-04	TP-1-EPH-2	2	1.16	9.7	8.78	89.2	
M2132-05	TP-1-VOC-TCL	3	1.13	9.71	8.76	88.9	
M2132-06	TP-2-COMP	4	1.15	9.97	8.78	86.5	
M2132-09	TP-2-EPH-2	5	1.17	9.67	7.92	79.4	
M2132-09D	TP-2-EPH-2DUP	6	1.14	9.69	7.94	79.5	
M2132-10	TP-2-VOC-TCL	7	1.16	9.98	9.11	90.1	
M2136-01	41763	8	1.13	9.69	8.49	86.0	
M2139-01	OIL-DEBRIS-1-2	54	1.00	2.00	2.00	100.0	OILY DEBRIS
M2140-01	203	9	1.12	9.51	8.57	88.8	
M2140-03	205	10	1.12	9.61	8.21	83.5	
M2141-01	202	11	1.13	9.86	8.62	85.8	
M2141-03	223	12	1.13	9.66	6.4	61.8	
M2141-05	935	13	1.14	9.82	8.85	88.8	
M2142-01	EO-04-042221	55	1.15	9.77	9.13	92.6	
M2142-02	EO-04-042221-EPH-2	56	1.12	9.7	9.1	93.0	
M2143-01	HD-01-042221	14	1.18	9.74	7.82	77.6	
M2143-02	HD-01-042221-EPH-2	15	1.18	9.86	8.77	87.4	
M2145-01	SU-02-042221	16	1.19	9.66	9.41	97.0	
M2145-02	SU-02-042221-EPH-2	17	1.13	9.57	9.47	98.8	
M2146-01	TR-04-042221	57	1.15	9.95	9.46	94.4	
M2146-02	TR-04-042221-EPH-2	58	1.18	9.69	8.79	89.4	
M2148-01	TP-18-17-COMP	18	1.18	9.74	7.39	72.5	
M2148-02	TP-18-17-EPH-2	19	1.15	9.95	7.42	71.3	
M2149-01	SB-36-0-2.0	20	1.14	9.96	8.84	87.3	
M2149-01D	SB-36-0-2.0DUP	21	1.13	9.94	8.86	87.7	
M2149-02	SB-36-COMP	22	1.13	9.83	8.99	90.3	
M2149-03	SB-36-10.5-11.0	23	1.12	9.63	7.83	78.8	



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 4/23/2021

OVENTEMP IN Celsius(°C): 106
Time IN: 17:00
In Date: 04/22/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:22
Out Date: 04/23/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB114113

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M2149-04	SB-37-0-2.0	24	1.1	9.86	9.04	90.6	
M2149-05	SB-37-COMP	25	1.14	9.92	8.91	88.5	
M2149-06	SB-37-9.5-10.0	26	1.16	9.92	7.67	74.3	
M2149-07	SB-38-0-2.0	27	1.12	9.58	7.82	79.2	
M2149-08	SB-38-COMP	28	1.13	9.72	7.94	79.3	
M2149-09	SB-38-9.5-10.0	29	1.18	9.55	7.73	78.3	
M2149-10	SB-39-0-2.0	30	1.19	9.64	8.66	88.4	
M2149-11	SB-39-COMP	31	1.15	9.6	8.26	84.1	
M2149-12	SB-39-9.5-10.0	32	1.15	9.97	8.28	80.8	
M2149-13	SB-40-0-2.0	33	1.14	9.94	9.17	91.3	
M2149-14	SB-40-COMP	34	1.13	9.57	8.52	87.6	
M2149-15	SB-40-9.5-10.0	35	1.14	9.94	8.68	85.7	
M2149-16	SB-40-9.5-10.0-DUP	36	1.17	9.97	8.38	81.9	
M2151-01	PE-1	37	1.11	9.51	8.45	87.4	
M2151-02	PE-2	38	1.15	9.66	8.45	85.8	
M2151-03	PE-3	39	1.16	9.62	8.39	85.5	
M2151-04	PE-4	40	1.14	9.98	9.27	92.0	
M2151-05	PE-5	41	1.16	9.54	8.52	87.8	
M2151-06	PE-6	42	1.1	9.61	8.56	87.7	
M2151-07	PE-7	43	1.11	9.8	8.49	84.9	
M2151-08	PE-8	44	1.14	9.84	9.26	93.3	
M2151-08D	PE-8DUP	45	1.19	9.86	9.28	93.3	
M2151-09	PE-9	46	1.19	9.78	9.18	93.0	
M2151-10	PE-10	47	1.1	9.71	8.87	90.2	
M2152-01	TP-3-COMP	48	1.14	9.94	8.29	81.3	
M2152-04	TP-3-EPH-2	49	1.13	9.63	8.1	82.0	
M2152-05	TP-3-VOC-TCL	50	1.14	9.82	7.2	69.8	
M2152-06	TP-4-COMP	51	1.18	9.84	8.41	83.5	



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 4/23/2021

OVENTEMP IN Celsius(°C): 106
Time IN: 17:00
In Date: 04/22/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:22
Out Date: 04/23/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB114113

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M2152-09	TP-4-EPH-2	52	1.18	9.75	8.52	85.6	
M2152-10	TP-4-VOC-TCL	53	1.19	9.91	7.97	77.8	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$