



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 4/26/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 17:20
In Date: 04/23/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:15
Out Date: 04/26/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB114134

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M2153-01	PLA-107	64	1.00	2.00	2.00	100.0	WIPE SAMPLE
M2153-02	PLA-108	65	1.00	2.00	2.00	100.0	WIPE SAMPLE
M2153-03	PLA-109	66	1.00	2.00	2.00	100.0	WIPE SAMPLE
M2153-04	PLA-110	67	1.00	2.00	2.00	100.0	WIPE SAMPLE
M2153-05	PLA-111	68	1.00	2.00	2.00	100.0	WIPE SAMPLE
M2153-06	PLA-112	69	1.00	2.00	2.00	100.0	WIPE SAMPLE
M2153-07	PLA-113	70	1.00	2.00	2.00	100.0	WIPE SAMPLE
M2153-08	PLA-114	71	1.00	2.00	2.00	100.0	WIPE SAMPLE
M2154-01	32722	72	1.00	2.00	2.00	100.0	OIL SAMPLE
M2155-01	25123	73	1.15	9.97	7.39	70.7	
M2157-01	KMA598M-END-1	1	1.00	2.00	2.00	100.0	PILC REEL
M2157-02	KMA598M-END-2	2	1.00	2.00	2.00	100.0	PILC REEL
M2157-03	HED302R-END-1	3	1.00	2.00	2.00	100.0	PILC REEL
M2157-04	HED302R-END-2	4	1.00	2.00	2.00	100.0	PILC REEL
M2158-01	KQA103H-END-1-1	5	1.00	2.00	2.00	100.0	PILC REEL
M2158-02	KQA103H-END-1-2	6	1.00	2.00	2.00	100.0	PILC REEL
M2158-03	KQA103H-END-1-3	7	1.00	2.00	2.00	100.0	PILC REEL
M2158-04	KQA103H-END-2-1	8	1.00	2.00	2.00	100.0	PILC REEL
M2158-05	KQA103H-END-2-2	9	1.00	2.00	2.00	100.0	PILC REEL
M2158-06	KQA103H-END-2-3	10	1.00	2.00	2.00	100.0	PILC REEL
M2158-07	KQA103H-END-3-1	11	1.00	2.00	2.00	100.0	PILC REEL
M2158-08	KQA103H-END-3-2	12	1.00	2.00	2.00	100.0	PILC REEL
M2158-09	KQA103H-END-3-3	13	1.00	2.00	2.00	100.0	PILC REEL
M2158-10	KQA103H-END-4-1	14	1.00	2.00	2.00	100.0	PILC REEL
M2158-11	KQA103H-END-4-2	15	1.00	2.00	2.00	100.0	PILC REEL
M2158-12	KQA103H-END-4-3	16	1.00	2.00	2.00	100.0	PILC REEL
M2159-01	BU-03-042321	74	1.16	9.82	9.19	92.7	
M2160-01	OK-03-042321	75	1.15	9.95	9.23	91.8	



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 4/26/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 17:20
In Date: 04/23/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:15
Out Date: 04/26/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB114134

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M2160-02	OK-03-042321-EPH-2	76	1.17	9.79	9.2	93.2	
M2163-01	TT-059-RW8-IDWSO-20210422-1	17	1.13	9.92	9.72	97.7	
M2163-02	TT-059-RW8-IDWSO-20210422-2	18	1.16	9.87	9.63	97.2	
M2163-03	TT-059-RW8-IDWSO-20210422-3	19	1.18	9.61	9.52	98.9	
M2164-01	BP-VPB-RW8-GW-243-245	20	1.15	9.97	5.32	47.3	SLUDGE SAMPLE
M2169-01	CF-1	21	1.14	9.61	8.27	84.2	
M2169-01D	CF-1DUP	22	1.18	9.63	8.29	84.1	
M2169-02	CF-2	23	1.19	9.82	8.53	85.1	
M2169-03	CF-3	24	1.19	9.84	8.64	86.1	
M2169-04	CF-4	25	1.19	9.74	8.37	84.0	
M2169-05	CF-5	26	1.16	9.74	8.6	86.7	
M2169-06	CF-6	27	1.13	9.76	8.54	85.9	
M2169-07	CF-7	28	1.19	9.55	8.53	87.8	
M2169-08	CF-8	29	1.1	9.56	8.36	85.8	
M2169-09	CF-9	30	1.19	9.91	8.45	83.3	
M2170-01	DS-1	31	1.13	9.82	8.99	90.4	
M2170-02	DN-1	32	1.16	9.51	8.64	89.6	
M2170-03	CS-1	33	1.19	9.71	8.53	86.2	
M2170-04	CN-1	34	1.19	9.67	8.04	80.8	
M2170-05	BN-1	35	1.17	9.6	8.65	88.7	
M2170-06	AN-1	36	1.15	9.97	9.26	92.0	
M2170-07	DS-1	37	1.13	9.52	8.73	90.6	
M2170-08	DS-1	38	1.13	9.52	8.73	90.6	
M2170-09	DN-1	39	1.19	9.63	8.97	92.2	
M2170-10	DN-1	40	1.19	9.63	8.97	92.2	
M2170-11	CS-1	41	1.14	9.94	8.86	87.7	
M2170-12	CS-1	42	1.14	9.94	8.86	87.7	



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 4/26/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 17:20
In Date: 04/23/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:15
Out Date: 04/26/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB114134

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M2170-13	CN-1	43	1.18	9.73	8.48	85.4	
M2170-14	CN-1	44	1.18	9.73	8.48	85.4	
M2170-15	BN-1	45	1.19	9.74	8.85	89.6	
M2170-16	BN-1	46	1.19	9.74	8.85	89.6	
M2170-17	AN-1	47	1.19	9.9	8.74	86.7	
M2170-18	AN-1	48	1.19	9.9	8.74	86.7	
M2171-01	P-1	49	1.17	9.72	8.22	82.5	
M2171-02	P-2	50	1.14	9.56	8.2	83.8	
M2171-03	T-1	51	1.1	9.74	8.34	83.8	
M2171-04	T-2	52	1.16	9.59	8.04	81.6	
M2171-05	T-3	53	1.19	9.74	8.2	82.0	
M2173-05	SVOC-GPC-BLANK	54	1.00	2.00	2.00	100.0	
M2173-06	PEST-GPC-BLANK	55	1.00	2.00	2.00	100.0	
M2173-07	PEST-GPC-BLANK-SPIKE	56	1.00	2.00	2.00	100.0	
M2173-08	PCB-GPC-BLANK	57	1.00	2.00	2.00	100.0	
M2173-09	PCB-GPC-BLANK-SPIKE	58	1.00	2.00	2.00	100.0	
M2173-10	SVOC-GCP2-BLANK	59	1.00	2.00	2.00	100.0	
M2173-11	PEST-GCP2-BLANK	60	1.00	2.00	2.00	100.0	
M2173-12	PEST-GCP2-BLANK-SPIKE	61	1.00	2.00	2.00	100.0	
M2173-13	PCB-GCP2-BLANK	62	1.00	2.00	2.00	100.0	
M2173-14	PCB-GCP2-BLANK-SPIKE	63	1.00	2.00	2.00	100.0	
M2176-01	STOCKPILE-VOC-1	77	1.14	9.74	8.86	89.8	
M2176-02	STOCKPILE-VOC-2	78	1.16	9.65	8.74	89.3	
M2176-03	STOCKPILE-VOC-3	79	1.16	9.92	8.91	88.5	
M2176-04	STOCKPILE-COMP-1	80	1.16	9.75	8.8	88.9	
M2176-06	STOCKPILE-VOC-4	81	1.12	9.64	8.87	91.0	
M2176-07	STOCKPILE-VOC-5	82	1.14	9.97	9.11	90.3	
M2176-08	STOCKPILE-VOC-6	83	1.16	9.88	9.13	91.4	



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 4/26/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 17:20
In Date: 04/23/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:15
Out Date: 04/26/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB114134

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M2176-09	STOCKPILE-COMP-2	84	1.16	9.71	8.81	89.5	
M2176-11	STOCKPILE-VOC-7	85	1.16	9.7	8.72	88.5	
M2176-12	STOCKPILE-VOC-8	86	1.15	9.94	9.14	90.9	
M2176-13	STOCKPILE-VOC-9	87	1.12	9.72	8.77	89.0	
M2176-14	STOCKPILE-COMP-3	88	1.14	9.97	9.14	90.6	
M2176-16	STOCKPILE-VOC-10	89	1.16	9.73	8.93	90.7	
M2176-17	STOCKPILE-VOC-11	90	1.1	9.83	9.06	91.2	
M2176-18	STOCKPILE-COMP-4	91	1.16	9.97	9.25	91.8	
M2179-01	N096253-02	92	1.15	9.98	1.68	6.0	SLUDGE SAMPLE
M2179-02	N096253-06	93	1.16	9.69	1.49	3.9	SLUDGE SAMPLE

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$