



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 5/25/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 16:30
In Date: 05/07/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 105
Time OUT: 10:00
Out Date: 05/08/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB114343

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
M2289-01	MH0AM0	1	1.14	9.94	7.44	71.6	
M2289-02	MH0AM0D	2	1.14	9.94	7.44	71.6	
M2289-03	MH0AM0S	3	1.14	9.94	7.44	71.6	
M2289-04	MH0AM1	4	1.16	9.74	7.27	71.2	
M2289-05	MH0AM2	5	1.19	9.66	6.84	66.7	
M2289-06	MH0AM3	6	1.17	9.98	6.66	62.3	
M2289-07	MH0AM4	7	1.12	9.63	5.83	55.3	
M2289-08	MH0AM5	8	1.18	9.97	4.96	43.0	
M2289-09	MH0AM6	9	1.13	9.52	6.91	68.9	
M2289-10	MH0AU9	10	1.16	9.97	7.04	66.7	
M2289-11	MH0AW0	11	1.16	9.91	7.94	77.5	
M2289-12	MH0AW1	12	1.1	9.89	8.07	79.3	
M2289-13	MH0AW2	13	1.16	9.97	7.89	76.4	
M2289-14	MH0AW3	14	1.16	9.98	8.29	80.8	
M2289-15	MH0AW4	15	1.12	9.76	8.33	83.4	
M2289-16	MH0AW5	16	1.15	9.88	6.63	62.8	
M2289-17	MH0AW6	17	1.15	9.97	6.92	65.4	
M2289-18	MH0AW7	18	1.16	9.92	7.7	74.7	
M2289-19	MH0AW8	19	1.15	9.97	7.09	67.3	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$