



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 5/11/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 16:30
In Date: 05/10/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:30
Out Date: 05/11/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB114370

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M2346-01	1-1000-GALLON-TANK	22	1.13	9.64	7.59	75.9	
M2346-02	2-1000-GALLON-TANK	23	1.14	9.72	7.22	70.9	
M2346-03	3-500-GALLON-TANK	24	1.16	9.75	8.68	87.5	
M2346-03D	3-500-GALLON-TANKDUP	25	1.17	9.77	8.7	87.6	
M2347-01	KMA7736-1-1	1	1.00	2.00	2.00	100.0	PILC REEL
M2347-02	KMA7736-1-2	2	1.00	2.00	2.00	100.0	PILC REEL
M2347-03	WOOD-REEL-1-1	3	1.00	2.00	2.00	100.0	PILC REEL
M2347-04	WOOD-REEL-1-2	4	1.00	2.00	2.00	100.0	PILC REEL
M2347-05	HEF090S-1-1	5	1.00	2.00	2.00	100.0	PILC REEL
M2347-06	HEF090S-1-2	6	1.00	2.00	2.00	100.0	PILC REEL
M2347-07	HEF090S-2-1	7	1.00	2.00	2.00	100.0	PILC REEL
M2347-08	HEF090S-2-2	8	1.00	2.00	2.00	100.0	PILC REEL
M2347-09	HEF090S-3-1	9	1.00	2.00	2.00	100.0	PILC REEL
M2347-10	HEF090S-3-2	10	1.00	2.00	2.00	100.0	PILC REEL
M2348-01	JEH669L-1-1	11	1.00	2.00	2.00	100.0	PILC REEL
M2348-02	JEH669L-1-2	12	1.00	2.00	2.00	100.0	PILC REEL
M2348-03	JEH669L-2-1	13	1.00	2.00	2.00	100.0	PILC REEL
M2348-04	JEH669L-2-2	14	1.00	2.00	2.00	100.0	PILC REEL
M2349-01	KQA110U-1-1	15	1.00	2.00	2.00	100.0	PILC REEL
M2349-02	KQA110U-1-2	16	1.00	2.00	2.00	100.0	PILC REEL
M2349-03	KQA110U-2-1	17	1.00	2.00	2.00	100.0	PILC REEL
M2349-04	KQA110U-2-2	18	1.00	2.00	2.00	100.0	PILC REEL
M2351-01	BC260909-1-1	19	1.00	2.00	2.00	100.0	PILC REEL
M2351-02	BC260909-1-2	20	1.00	2.00	2.00	100.0	PILC REEL
M2351-03	BC260909-1-3	21	1.00	2.00	2.00	100.0	PILC REEL

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$