



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 5/12/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 05/11/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:26
Out Date: 05/12/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB114385

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M2354-01	EO-01-051121	1	1.18	9.68	8.37	84.6	
M2354-02	EO-01-051121-EPH-2	2	1.14	9.9	8.65	85.7	
M2355-01	JEH669L-END-1-1	3	1.00	2.00	2.00	100.0	PILC REEL
M2355-02	JEH669L-END-1-2	4	1.00	2.00	2.00	100.0	PILC REEL
M2355-03	KZH012P-END-1-1	5	1.00	2.00	2.00	100.0	PILC REEL
M2355-04	KZH012P-END-1-2	6	1.00	2.00	2.00	100.0	PILC REEL
M2356-01	OR-03-051121	7	1.18	9.67	8.91	91.0	
M2356-02	OR-03-051121-EPH-2	8	1.19	9.52	8.79	91.2	
M2361-01	CL-02-051121	9	1.18	9.83	8.6	85.8	
M2361-02	CL-02-051121-EPH-2	10	1.12	9.87	9.00	90.1	
M2362-01	149TH-ST-DRUM-1-PPE	11	1.00	2.00	2.00	100.0	plastic sample
M2362-02	TREMONT-AVE-DRUM-8-PPE	12	1.00	2.00	2.00	100.0	plastic sample
M2362-03	TREMONT-AVE-DRUM-9-POLY	13	1.00	2.00	2.00	100.0	plastic sample
M2364-01	AST2-COPR4-1	14	1.18	9.71	8.62	87.2	
M2364-01D	AST2-COPR4-1DUP	15	1.18	9.73	8.64	87.3	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$