



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 5/13/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 16:42
In Date: 05/12/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:22
Out Date: 05/13/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB114407

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M2366-01	JG-SPOILS-5-11-21-A	1	1.15	9.96	6.37	59.3	
M2366-02	JG-SPOILS-5-11-21-B	2	1.1	9.87	6.21	58.3	
M2366-03	JG-SPOILS-5-11-21-C	3	1.18	9.85	5.76	52.8	
M2366-03D	JG-SPOILS-5-11-21-CDUP	4	1.11	9.87	5.78	53.3	
M2367-01	PV-SP1	23	1.15	9.97	9.18	91.0	
M2367-03	PV-SP1-EPH-2	24	1.15	9.97	9.38	93.3	
M2367-04	PV-SP1-EPH-3	25	1.16	9.94	9.28	92.5	
M2367-05	PV-SP2	26	1.12	9.91	8.77	87.0	
M2367-07	PV-SP2-EPH-2	27	1.15	9.97	8.74	86.1	
M2367-08	PV-SP2-EPH-3	28	1.16	9.74	8.52	85.8	
M2369-01	RO-365	5	1.15	9.68	8.92	91.1	
M2369-02	RO-365	6	1.16	9.54	8.98	93.3	
M2370-01	367	29	1.19	9.62	8.72	89.3	
M2370-02	367	30	1.19	9.71	8.54	86.3	
M2371-01	VNJ-206	7	1.15	9.97	8.87	87.5	
M2372-01	CHRT-24702	8	1.16	9.81	8.79	88.2	
M2372-02	CHRT-24702-EPH-2	9	1.14	9.61	8.84	90.9	
M2372-03	CHRT-27051	10	1.16	9.54	8.93	92.7	
M2372-04	CHRT-27051-EPH-2	11	1.15	9.96	8.89	87.9	
M2372-05	CHRT-26598	12	1.19	9.87	9.42	94.8	
M2372-06	CHRT-26598-EPH-2	13	1.16	9.51	9.33	97.8	
M2372-07	RBR-200032	14	1.16	9.84	7.94	78.1	
M2372-08	RBR-200032-EPH-2	15	1.17	9.61	7.29	72.5	
M2372-09	CHRT-26516	16	1.1	9.62	8.58	87.8	
M2372-10	CHRT-26516-EPH-2	17	1.19	9.81	8.71	87.2	
M2372-11	CHRT-26855	18	1.15	9.94	7.06	67.2	
M2372-12	CHRT-26855-EPH-2	19	1.16	9.9	7.92	77.3	
M2372-13	TOTE-3	20	1.00	2.00	2.00	100.0	OIL SAMPLE



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 5/13/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 16:42
In Date: 05/12/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:22
Out Date: 05/13/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB114407

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M2373-01	OR-02-051221	31	1.15	9.76	9.37	95.5	
M2373-01D	OR-02-051221DUP	32	1.19	9.78	9.39	95.5	
M2373-02	OR-02-051221-EPH-2	33	1.17	9.68	9.23	94.7	
M2375-01	SU-03-051221	21	1.16	9.81	9.34	94.6	
M2375-02	SU-03-051221-EPH-2	22	1.17	9.79	9.39	95.4	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$