



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 5/24/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 05/13/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:26
Out Date: 05/14/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB114436

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M2383-01	KSB-01 (6-8)	1	1.15	9.61	8.47	86.5	
M2383-02	KSB-02 (7.5-9.5)	2	1.13	9.82	7.4	72.2	
M2383-03	KSB-03 (7-9)	3	1.16	9.71	8.79	89.2	
M2383-04	M2383-03MS	4	1.16	9.71	8.79	89.2	
M2383-05	M2383-03MSD	5	1.16	9.71	8.79	89.2	
M2383-06	DUP-01-05112021	6	1.13	9.66	7.15	70.6	
M2385-01	33077	10	1.00	2.00	2.00	100.0	gasoline
M2385-02	32798	11	1.00	2.00	2.00	100.0	OILY DEBRIS
M2386-01	32519	12	1.12	9.56	8.65	89.2	
M2387-01	41789	15	1.12	9.54	9.49	99.4	
M2387-02	30782	16	1.00	2.00	2.00	100.0	OILY DEBRIS
M2388-01	41785-LIGHT	17	1.00	2.00	2.00	100.0	OIL SAMPLE
M2388-03	41785-DARK	18	1.00	2.00	2.00	100.0	OIL SAMPLE
M2388-05	41786-LIGHT	19	1.00	2.00	2.00	100.0	OIL SAMPLE
M2388-07	41786-DARK	20	1.00	2.00	2.00	100.0	OIL SAMPLE
M2389-01	BURNETT-AVE-1-3	41	1.15	9.95	9.34	93.1	
M2389-02	BURNETT-AVE-1-3	42	1.16	9.98	9.11	90.1	
M2389-03	BURNETT-AVE-4-6	43	1.15	9.55	8.65	89.3	
M2389-04	BURNETT-AVE-4-6	44	1.16	9.62	8.75	89.7	
M2389-05	CASINO-AVE-1-2	45	1.19	9.57	8.64	88.9	
M2389-06	CASINO-AVE-1-2	46	1.12	9.81	8.69	87.1	
M2389-07	CASINO-AVE-4-6	47	1.12	9.92	8.17	80.1	
M2389-08	CASINO-AVE-4-6	48	1.16	9.79	8.2	81.6	
M2389-08D	CASINO-AVE-4-6DUP	49	1.17	9.81	8.22	81.6	
M2390-01	01-A-01-B-01-C	7	1.00	2.00	2.00	100.0	CAULKING SAMPLE
M2390-02	02-A-02-B-02-C	8	1.00	2.00	2.00	100.0	CAULKING SAMPLE
M2390-03	03-A-03-B-03-C	9	1.00	2.00	2.00	100.0	CAULKING SAMPLE
M2391-01	HEF239Y-2-1	21	1.00	2.00	2.00	100.0	PILC REEL



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 5/24/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 05/13/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:26
Out Date: 05/14/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB114436

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
M2391-02	HEF239Y-2-2	22	1.00	2.00	2.00	100.0	PILC REEL
M2391-03	HEF239Y-3-1	23	1.00	2.00	2.00	100.0	PILC REEL
M2391-04	HEF239Y-3-2	24	1.00	2.00	2.00	100.0	PILC REEL
M2392-01	SB-01-3.5-4.0	25	1.16	9.74	8.73	88.2	
M2392-02	SB-01-COMP	26	1.13	9.78	8.74	88.0	
M2392-04	SB-02-3.5-4.0	27	1.16	9.72	8.4	84.6	
M2392-05	SB-02-COMP	28	1.12	9.97	8.73	86.0	
M2392-07	SB-03-3.5-4.0	29	1.14	9.96	8.76	86.4	
M2392-08	SB-03-COMP	30	1.16	9.61	8.35	85.1	
M2392-10	SB-04-4.5-5.0	31	1.18	9.95	8.8	86.9	
M2392-11	SB-04-COMP	32	1.19	9.88	8.96	89.4	
M2392-13	SB-05-4.0-4.5	33	1.15	9.79	7.96	78.8	
M2392-14	SB-05-COMP	34	1.15	9.67	7.87	78.9	
M2392-16	SB-06-4.5-5.0	35	1.17	9.89	8.62	85.4	
M2392-17	SB-06-COMP	36	1.15	9.96	8.53	83.8	
M2392-19	SB-07-4.5-5.0	37	1.16	9.72	8.81	89.4	
M2392-20	SB-07-COMP	38	1.13	9.88	8.99	89.8	
M2392-22	SB-08-3.5-4.0	39	1.18	9.66	8.81	90.0	
M2392-23	SB-08-COMP	40	1.12	9.56	8.78	90.8	
M2403-01	PE-5- (5-5.5)	13	1.14	9.54	7.92	80.7	
M2403-01D	PE-5- (5-5.5) DUP	14	1.16	9.56	7.95	80.8	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$