



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 5/14/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 14:00
In Date: 05/13/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:00
Out Date: 05/14/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB114439

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M2400-13	X3373	1	1.00	2.00	2.00	100.0	P.E. SAMPLE
M2400-14	X3373D	2	1.00	2.00	2.00	100.0	P.E. SAMPLE
M2400-15	X3373S	3	1.00	2.00	2.00	100.0	P.E. SAMPLE
M2400-16	X3374	4	1.00	2.00	2.00	100.0	P.E. SAMPLE
M2400-17	X3374D	5	1.00	2.00	2.00	100.0	P.E. SAMPLE
M2400-18	X3374S	6	1.00	2.00	2.00	100.0	P.E. SAMPLE
M2400-19	X3375	7	1.00	2.00	2.00	100.0	P.E. SAMPLE
M2400-20	X3375D	8	1.00	2.00	2.00	100.0	P.E. SAMPLE
M2400-21	X3375S	9	1.00	2.00	2.00	100.0	P.E. SAMPLE
M2400-22	X3376	10	1.00	2.00	2.00	100.0	P.E. SAMPLE
M2400-23	X3376D	11	1.00	2.00	2.00	100.0	P.E. SAMPLE
M2400-24	X3376S	12	1.00	2.00	2.00	100.0	P.E. SAMPLE

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$