



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 5/18/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 16:30
In Date: 05/17/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:00
Out Date: 05/18/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB114493

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M2405-01	BUR-1013	1	1.00	2.00	2.00	100.0	OILY DEBRIS
M2405-02	BUR-1019	2	1.00	2.00	2.00	100.0	OILY DEBRIS
M2405-03	BUR-1013-1019	3	1.00	2.00	2.00	100.0	OILY DEBRIS
M2420-01	CONCESSION-CP-1	4	1.15	9.94	9.24	92.0	
M2420-02	CONCESSION-CP-2	5	1.12	9.67	9.37	96.5	
M2420-03	CONCESSION-CP-3	6	1.13	9.84	9.4	94.9	
M2420-04	CONCESSION-CP-4	7	1.13	9.97	9.44	94.0	
M2420-05	CONCESSION-CP-5	8	1.16	9.87	9.41	94.7	
M2420-06	PLAYGROUND-CP-1	9	1.15	9.76	4.59	40.0	
M2421-01	K328-1A	10	1.13	9.59	8.67	89.1	
M2421-02	K328-DUP1	11	1.15	9.5	8.4	86.8	
M2421-04	K328-2A	12	1.16	9.94	9.27	92.4	
M2421-06	K328-3A	13	1.13	9.61	8.63	88.4	
M2421-08	K328-4A	14	1.16	9.57	8.97	92.9	
M2421-10	K328-5A	15	1.18	9.76	8.76	88.3	
M2421-12	K328-6A	16	1.17	9.62	9.03	93.0	
M2421-13	K328-DUP2	17	1.12	9.6	8.95	92.3	
M2421-15	K328-7A	18	1.14	9.67	8.59	87.3	
M2421-17	K328-8A	19	1.15	9.88	8.84	88.1	
M2421-19	K328-9A	20	1.16	9.79	8.85	89.1	
M2422-01	K328-10A	21	1.14	9.82	8.99	90.4	
M2422-01D	K328-10ADUP	22	1.15	9.85	9.02	90.5	
M2422-03	K328-11A	23	1.17	9.6	9.00	92.9	
M2422-05	K328-DUP3	24	1.15	9.96	9.44	94.1	
M2422-06	K328-12A	25	1.16	9.78	8.02	79.6	
M2422-09	K328-13A	26	1.15	9.96	8.56	84.1	
M2422-12	K328-14A	27	1.19	9.64	8.34	84.6	
M2422-15	K328-15A	28	1.16	9.68	8.06	81.0	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 5/18/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 16:30
In Date: 05/17/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:00
Out Date: 05/18/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB114493

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
M2422-18	K328-DUP4	29	1.16	9.67	8.25	83.3	
M2422-19	K328-16A	30	1.13	9.58	7.84	79.4	
M2423-02	K328-17A	31	1.13	9.77	8.05	80.1	
M2423-03	K328-DUP5	32	1.13	9.84	8.06	79.6	
M2423-06	K328-18A	33	1.1	9.71	3.99	33.6	
M2423-07	K328-19A	34	1.14	9.67	4.63	40.9	
M2423-08	K328-20A	35	1.12	9.87	4.74	41.4	
M2423-09	K328-21A	36	1.15	9.56	4.43	39.0	
M2423-10	K328-22A	37	1.16	9.82	3.92	31.9	
M2423-11	K328-23A	38	1.14	9.94	4.14	34.1	
M2423-12	K328-24A	39	1.18	9.86	6.16	57.4	
M2423-13	K328-25A	40	1.13	9.57	5.24	48.7	
M2423-14	K328-26A	41	1.13	9.82	5.05	45.1	
M2423-15	K328-27A	42	1.16	9.93	5.22	46.3	
M2423-16	K328-28A	43	1.16	9.51	4.27	37.2	
M2423-17	K328-DUP6	44	1.16	9.74	4.67	40.9	
M2423-18	K328-29A	45	1.13	9.6	4.9	44.5	
M2423-19	K328-30A	46	1.16	9.54	4.15	35.7	
M2423-20	K328-31A	47	1.1	9.78	5.08	45.9	
M2423-21	K328-32A	48	1.11	9.56	4.19	36.4	
M2423-22	K328-33A	49	1.16	9.5	4.4	38.8	
M2423-23	K328-34A	50	1.16	9.74	4.00	33.1	
M2423-23D	K328-34ADUP	51	1.15	9.76	4.15	34.8	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$