



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 6/8/2021

OVENTEMP IN Celsius(°C): 108
Time IN: 17:00
In Date: 06/07/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:13
Out Date: 06/08/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB114824

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M2625-01	18-B12 (138-142)	1	1.14	9.74	8.03	80.1	
M2625-02	18-B12 (142-146)	2	1.12	9.88	8.02	78.8	
M2625-03	18-B12 (146-150)	3	1.14	9.63	7.98	80.6	
M2625-04	18-B12 (150-154)	4	1.15	9.97	8.46	82.9	
M2625-05	18-B12 (154-158)	5	1.16	9.54	6.54	64.2	
M2625-06	18-B12 (158-162)	6	1.17	9.84	8.43	83.7	
M2625-07	18-B12 (162-166)	7	1.14	9.86	8.27	81.8	
M2625-08	18-B12 (166-170)	8	1.12	9.77	7.93	78.7	
M2625-09	18-B12 (170-174)	9	1.15	9.54	8.23	84.4	
M2625-10	18-B12 (174-178)	10	1.1	9.88	8.64	85.9	
M2625-11	18-B12 (178-182)	11	1.16	9.97	8.88	87.6	
M2625-12	18-B12 (182-186)	12	1.14	9.95	9.00	89.2	
M2625-13	18-B12 (186-190)	13	1.14	9.91	8.47	83.6	
M2625-14	18-B12 (190-194)	14	1.18	9.7	8.32	83.8	
M2625-15	18-B12 (194-198)	15	1.16	9.66	8.44	85.6	
M2625-16	18-B12 (198-202)	16	1.16	9.97	8.64	84.9	
M2625-17	18-B12 (202-206)	17	1.12	9.62	8.54	87.3	
M2625-18	18-B12 (206-210)	18	1.15	9.97	9.3	92.4	
M2625-19	18-B12 (210-214)	19	1.13	9.55	8.23	84.3	
M2625-20	18-B12 (214-218)	20	1.13	9.59	8.2	83.6	
M2625-20D	18-B12 (214-218) DUP	21	1.15	9.6	8.22	83.7	
M2626-01	18-B12 (218-222)	22	1.14	9.97	9.04	89.5	
M2626-02	18-B12 (222-226)	23	1.15	9.95	7.89	76.6	
M2626-03	18-B12 (226-230)	24	1.15	9.66	7.71	77.1	
M2626-04	18-B12-DUP-01	25	1.15	9.97	8.49	83.2	
M2626-05	18-B12 (X)	26	1.13	9.91	8.72	86.4	
M2626-06	18-B12 (XX)	27	1.13	9.7	8.33	84.0	
M2628-01	HD-02-060721	31	1.12	9.9	9.06	90.4	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 6/8/2021

OVENTEMP IN Celsius(°C): 108
Time IN: 17:00
In Date: 06/07/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:13
Out Date: 06/08/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB114824

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M2628-02	HD-02-060721-EPH-2	32	1.13	9.76	8.96	90.7	
M2629-01	OK-02-060721	29	1.14	9.71	9.47	97.2	
M2629-02	OK-02-060721-EPH-2	30	1.13	9.85	9.55	96.6	
M2633-01	WASTE	28	1.16	9.71	9.52	97.8	
M2634-01	R44-STONE-1	33	1.00	2.00	2.00	100.0	STONE SAMPLE,100% SOLIDS
M2634-02	R44-STONE-1-EPH-2	34	1.00	2.00	2.00	100.0	STONE SAMPLE,100% SOLIDS
M2634-04	R44-STONE-2	35	1.00	2.00	2.00	100.0	STONE SAMPLE,100% SOLIDS
M2634-05	R44-STONE-2-EPH-2	36	1.00	2.00	2.00	100.0	STONE SAMPLE,100% SOLIDS
M2637-01	18-B12 (230-234)	37	1.14	9.98	7.91	76.6	
M2637-02	18-B12 (234-238)	38	1.16	9.96	8.72	85.9	
M2637-03	18-B12 (238-242)	39	1.12	9.69	7.84	78.4	
M2637-04	18-B12 (242-246)	40	1.13	9.71	8.29	83.4	
M2637-05	18-B12 (246-248)	41	1.1	9.86	8.39	83.2	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$