



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 6/16/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 16:40
In Date: 06/15/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:10
Out Date: 06/16/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB114965

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M2711-01	PE-1	1	1.15	9.98	9.00	88.9	
M2711-02	PE-2	2	1.15	9.95	8.83	87.3	
M2711-03	PE-3	3	1.12	9.63	8.93	91.8	
M2711-04	PE-4	4	1.16	9.62	7.91	79.8	
M2711-05	PE-B	5	1.13	9.97	8.73	86.0	
M2714-01	33241	41	1.00	2.00	2.00	100.0	DEBRIS
M2715-01	33248	42	1.15	9.71	9.7	99.9	
M2715-03	33258	43	1.00	2.00	2.00	100.0	gasoline
M2716-01	33254	44	1.11	9.72	9.2	94.0	
M2718-01	WA-1	12	1.00	2.00	2.00	100.0	WIPE SAMPLE,100% SOLIDS
M2718-02	WA-2	13	1.00	2.00	2.00	100.0	WIPE SAMPLE,100% SOLIDS
M2718-03	WA-3	14	1.00	2.00	2.00	100.0	WIPE SAMPLE,100% SOLIDS
M2718-04	WA-4	15	1.00	2.00	2.00	100.0	WIPE SAMPLE,100% SOLIDS
M2718-05	WA-5	16	1.00	2.00	2.00	100.0	WIPE SAMPLE,100% SOLIDS
M2718-06	WA-6	17	1.00	2.00	2.00	100.0	WIPE SAMPLE,100% SOLIDS
M2718-07	WA-7	18	1.00	2.00	2.00	100.0	WIPE SAMPLE,100% SOLIDS
M2718-08	WA-8	19	1.00	2.00	2.00	100.0	WIPE SAMPLE,100% SOLIDS
M2718-09	WA-9	20	1.00	2.00	2.00	100.0	WIPE SAMPLE,100% SOLIDS
M2718-10	WA-10	21	1.00	2.00	2.00	100.0	WIPE SAMPLE,100% SOLIDS
M2718-11	WA-11	22	1.00	2.00	2.00	100.0	WIPE SAMPLE,100% SOLIDS
M2718-12	WA-12	23	1.00	2.00	2.00	100.0	WIPE SAMPLE,100% SOLIDS
M2718-13	WA-13	24	1.00	2.00	2.00	100.0	WIPE SAMPLE,100% SOLIDS
M2718-14	WA-14	25	1.00	2.00	2.00	100.0	WIPE SAMPLE,100% SOLIDS
M2718-15	WA-15	26	1.00	2.00	2.00	100.0	WIPE SAMPLE,100% SOLIDS
M2718-16	WA-16	27	1.00	2.00	2.00	100.0	WIPE SAMPLE,100% SOLIDS
M2718-17	WA-17	28	1.00	2.00	2.00	100.0	WIPE SAMPLE,100% SOLIDS
M2718-18	WA-18	29	1.00	2.00	2.00	100.0	WIPE SAMPLE,100% SOLIDS
M2718-19	WA-19	30	1.00	2.00	2.00	100.0	WIPE SAMPLE,100% SOLIDS



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 6/16/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 16:40
In Date: 06/15/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:10
Out Date: 06/16/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB114965

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
M2718-20	WA-20	31	1.00	2.00	2.00	100.0	WIPE SAMPLE,100% SOLIDS
M2718-21	WA-21	32	1.00	2.00	2.00	100.0	WIPE SAMPLE,100% SOLIDS
M2718-22	WA-22	33	1.00	2.00	2.00	100.0	WIPE SAMPLE,100% SOLIDS
M2718-23	WA-23	34	1.00	2.00	2.00	100.0	WIPE SAMPLE,100% SOLIDS
M2718-24	WA-24	35	1.00	2.00	2.00	100.0	WIPE SAMPLE,100% SOLIDS
M2718-25	WA-25	36	1.00	2.00	2.00	100.0	WIPE SAMPLE,100% SOLIDS
M2718-26	WA-26	37	1.00	2.00	2.00	100.0	WIPE SAMPLE,100% SOLIDS
M2718-27	WA-27	38	1.00	2.00	2.00	100.0	WIPE SAMPLE,100% SOLIDS
M2718-28	WA-28	39	1.00	2.00	2.00	100.0	WIPE SAMPLE,100% SOLIDS
M2718-29	WA-29	40	1.00	2.00	2.00	100.0	WIPE SAMPLE,100% SOLIDS
M2719-01	OK-03-061521	6	1.12	9.55	9.22	96.1	
M2719-02	OK-03-061521-EPH-2	7	1.13	9.89	9.57	96.3	
M2722-01	P001-SS181-0001-01	8	1.16	9.74	8.96	90.9	
M2722-02	P001-SS181-0001-02	9	1.19	9.72	8.83	89.6	
M2722-03	P001-SS182-0001-01	10	1.12	9.66	8.81	90.0	
M2722-03D	P001-SS182-0001-01DUP	11	1.12	9.68	8.85	90.3	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$