



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 6/18/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 16:30
In Date: 06/17/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:13
Out Date: 06/18/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB114995

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M2747-01	A1 (WINDROW-LEFT-OF-ENTRANCE)	1	1.14	9.96	9.25	92.0	
M2747-02	A2 (WINDROW-MIDDLE)	2	1.16	9.75	9.09	92.3	
M2747-03	A3 (WINDROW-RIGHT)	3	1.14	9.94	8.76	86.6	
M2747-04	X4 (TOPSOIL-PILE)	4	1.16	9.83	9.09	91.5	
M2747-05	B1 (WINDROW-1-ABOVE)	5	1.15	9.95	9.37	93.4	
M2747-06	B2 (WINDROW-2-ABOVE)	6	1.18	9.91	9.47	95.0	
M2747-07	B3 (WINDROW-3-ABOVE)	7	1.12	9.76	9.24	94.0	
M2747-08	X5 (TOPSOIL-PILE)	8	1.15	9.97	9.45	94.1	
M2747-09	C1 (WINDROW-LEFT-OF-ENTRANCE)	9	1.16	9.63	9.35	96.7	
M2747-10	C2 (WINDROW-MIDDLE)	10	1.14	9.98	9.39	93.3	
M2747-11	C33 (WINDROW-RIGHT)	11	1.12	9.95	9.38	93.5	
M2747-12	X6 (TOPSOIL-PILE)	12	1.12	9.62	8.00	80.9	
M2747-12D	X6 (TOPSOIL-PILE) DUP	13	1.18	9.64	8.02	80.9	
M2749-01	PX-01_14.5-15.0	14	1.14	9.98	8.24	80.3	
M2749-02	PX-02_14.5-15.0	15	1.15	9.95	7.5	72.2	
M2749-03	PX-03_14.5-15.0	16	1.11	9.77	8.04	80.0	
M2749-04	PX-04_14.5-15.0	17	1.13	9.77	8.53	85.6	
M2749-05	PP-01_4.5-5.0	18	1.15	9.74	8.2	82.1	
M2749-06	RF-01_3.5-4.0	19	1.17	9.98	8.11	78.8	
M2754-01	41885	49	1.00	2.00	2.00	100.0	WIPE SAMPLE
M2754-02	41886	50	1.00	2.00	2.00	100.0	WIPE SAMPLE
M2754-03	41887	51	1.00	2.00	2.00	100.0	WIPE SAMPLE
M2755-01	HEE262Y-N-1	25	1.00	2.00	2.00	100.0	PILC REEL
M2755-02	HEE262Y-N-2	26	1.00	2.00	2.00	100.0	PILC REEL
M2756-01	230215-N-1-1	27	1.00	2.00	2.00	100.0	PILC REEL
M2756-02	230215-N-1-2	28	1.00	2.00	2.00	100.0	PILC REEL
M2756-03	230215-N-2-1	29	1.00	2.00	2.00	100.0	PILC REEL



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 6/18/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 16:30
In Date: 06/17/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:13
Out Date: 06/18/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB114995

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M2756-04	230215-N-2-2	30	1.00	2.00	2.00	100.0	PILC REEL
M2756-05	81644-N-1	31	1.00	2.00	2.00	100.0	PILC REEL
M2756-06	81644-N-2	32	1.00	2.00	2.00	100.0	PILC REEL
M2756-07	HEF317G-N-1	33	1.00	2.00	2.00	100.0	PILC REEL
M2756-08	HEF317G-N-2	34	1.00	2.00	2.00	100.0	PILC REEL
M2756-09	HEE9091-N-1	35	1.00	2.00	2.00	100.0	PILC REEL
M2756-10	HEE9091-N-2	36	1.00	2.00	2.00	100.0	PILC REEL
M2757-01	KZA467N-N-1	37	1.00	2.00	2.00	100.0	PILC REEL
M2757-02	KZA467N-N-2	38	1.00	2.00	2.00	100.0	PILC REEL
M2757-03	JEB542P-N-1	39	1.00	2.00	2.00	100.0	PILC REEL
M2757-04	JEB542P-N-2	40	1.00	2.00	2.00	100.0	PILC REEL
M2757-05	JEH851S-N-1	41	1.00	2.00	2.00	100.0	PILC REEL
M2757-06	JEH851S-N-2	42	1.00	2.00	2.00	100.0	PILC REEL
M2757-07	JEH705N-N-1	43	1.00	2.00	2.00	100.0	PILC REEL
M2757-08	JEH705N-N-2	44	1.00	2.00	2.00	100.0	PILC REEL
M2757-09	FDA749E-N-1	45	1.00	2.00	2.00	100.0	PILC REEL
M2757-10	FDA749E-N-2	46	1.00	2.00	2.00	100.0	PILC REEL
M2757-11	HED184J-N-1	47	1.00	2.00	2.00	100.0	PILC REEL
M2757-12	HED184J-N-2	48	1.00	2.00	2.00	100.0	PILC REEL
M2759-01	TR-05-061721	52	1.14	9.7	8.47	85.6	
M2759-02	TR-05-061721-EPH-2	53	1.14	9.94	8.92	88.4	
M2760-01	01-A-01-B-01-C	21	1.00	2.00	2.00	100.0	CAULKING SAMPLE
M2760-02	02-A-02-B-02-C	22	1.00	2.00	2.00	100.0	CAULKING SAMPLE
M2762-01	WC_SOLIDS_DBD	20	1.18	9.72	9.33	95.4	
M2764-01	PHC-250EM-001	23	1.11	9.94	9.37	93.5	
M2764-03	PHC-250EM-002	24	1.12	9.79	8.59	86.2	
M2768-01	OK-01-061721	54	1.18	9.67	9.28	95.4	
M2768-02	OK-01-061721-EPH-2	55	1.15	9.79	9.46	96.2	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$