



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 6/28/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 06/25/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:11
Out Date: 06/26/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB115106

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M2856-01	SB-01-9.5-10.0	1	1.14	9.68	9.27	95.2	
M2856-02	SB-01-COMP	2	1.12	9.97	9.3	92.4	
M2856-03	SB-02-9.5-10.0	3	1.15	9.59	8.83	91.0	
M2856-04	SB-02-COMP	4	1.13	9.8	9.22	93.3	
M2856-05	SB-03-9.5-10.0	5	1.16	9.82	9.34	94.5	
M2856-06	SB-03-COMP	6	1.16	9.87	9.16	91.8	
M2856-07	SB-04-9.5-10.0	7	1.18	9.65	9.26	95.4	
M2856-08	SB-04-COMP	8	1.16	9.88	9.44	95.0	
M2856-09	SB-05-9.5-10.0	9	1.15	9.6	9.21	95.4	
M2856-10	SB-05-COMP	10	1.12	9.93	9.47	94.8	
M2856-11	SB-06-9.5-10.0	11	1.14	9.94	9.28	92.5	
M2856-12	SB-06-COMP	12	1.12	9.62	9.17	94.7	
M2856-13	SB-07-9.5-10.0	13	1.15	9.97	9.1	90.1	
M2856-14	SB-07-COMP	14	1.16	9.59	9.24	95.8	
M2856-15	SB-08-9.5-10.0	15	1.15	9.66	8.19	82.7	
M2856-16	SB-08-COMP	16	1.15	9.94	8.39	82.4	
M2856-17	SB-09-9.5-10.0	17	1.13	9.81	9.44	95.7	
M2856-18	SB-09-COMP	18	1.14	9.8	9.06	91.5	
M2856-19	SB-10-9.5-10.0	19	1.12	9.69	9.09	93.0	
M2856-20	SB-10-COMP	20	1.17	9.51	9.34	98.0	
M2857-01	SB-11-9.5-10.0	21	1.16	9.93	9.44	94.4	
M2857-02	SB-11-COMP	22	1.1	9.82	9.14	92.2	
M2857-03	SB-12-9.5-10.0	23	1.15	9.96	9.48	94.6	
M2857-04	SB-12-COMP	24	1.16	9.97	8.85	87.3	
M2857-05	SB-13-9.5-10.0	25	1.13	9.51	8.59	89.0	
M2857-06	SB-13-COMP	26	1.14	9.95	9.34	93.1	
M2857-07	SB-14-9.5-10.0	27	1.11	9.89	9.48	95.3	
M2857-08	SB-14-COMP	28	1.15	9.98	9.1	90.0	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 6/28/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 06/25/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:11
Out Date: 06/26/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB115106

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M2857-09	SB-15-9.5-10.0	29	1.14	9.92	8.65	85.5	
M2857-10	SB-15-COMP	30	1.15	9.52	8.19	84.1	
M2857-11	SB-16-9.5-10.0	31	1.15	9.9	9.06	90.4	
M2857-12	SB-16-COMP	32	1.15	9.96	9.28	92.3	
M2857-13	SB-17-9.5-10.0	33	1.15	9.97	9.34	92.9	
M2857-14	SB-17-COMP	34	1.14	9.97	9.48	94.5	
M2857-15	SB-18-9.5-10.0	35	1.14	9.95	9.47	94.6	
M2857-16	SB-18-COMP	36	1.13	9.62	7.66	76.9	
M2859-01	PILE-1-062421	37	1.18	9.56	9.41	98.2	
M2859-02	PILE-2-062421	38	1.16	9.87	9.31	93.6	
M2864-01	VNJ-206	56	1.18	9.92	9.1	90.6	
M2864-02	VNJ-206	57	1.15	9.74	8.42	84.6	
M2864-03	12002	58	1.18	9.7	8.55	86.5	
M2864-04	12002	59	1.13	9.93	8.03	78.4	
M2865-01	OR-01-062521	60	1.15	9.78	9.34	94.9	
M2865-02	OR-01-062521-EPH-2	61	1.15	9.58	9.28	96.4	
M2865-02D	OR-01-062521-EPH-2DUP	62	1.14	9.64	9.3	96.0	
M2866-01	124055	49	1.13	9.52	9.07	94.6	
M2866-03	124030	50	1.11	9.65	8.32	84.4	
M2866-05	124019	51	1.12	9.69	8.47	85.8	
M2866-07	124042	52	1.15	9.97	9.14	90.6	
M2867-01	OILY-SOCK-DRUM	53	1.00	2.00	2.00	100.0	oily-stock-drum
M2867-03	OILY-STONE-DRUM	54	1.00	2.00	2.00	100.0	oily-stone-drum
M2868-01	OILY-DEBRIS-TRE-1169	55	1.00	2.00	2.00	100.0	OILY DEBRIS
M2872-05	SVOC-GPC-BLANK	39	1.00	2.00	2.00	100.0	
M2872-06	PEST-GPC-BLANK	40	1.00	2.00	2.00	100.0	
M2872-07	PEST-GPC-BLANK-SPIKE	41	1.00	2.00	2.00	100.0	
M2872-08	PCB-GPC-BLANK	42	1.00	2.00	2.00	100.0	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 6/28/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 06/25/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:11
Out Date: 06/26/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB115106

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M2872-09	PCB-GPC-BLANK-SPIKE	43	1.00	2.00	2.00	100.0	
M2872-10	SVOC-GCP2-BLANK	44	1.00	2.00	2.00	100.0	
M2872-11	PEST-GCP2-BLANK	45	1.00	2.00	2.00	100.0	
M2872-12	PEST-GCP2-BLANK-SPIKE	46	1.00	2.00	2.00	100.0	
M2872-13	PCB-GCP2-BLANK	47	1.00	2.00	2.00	100.0	
M2872-14	PCB-GCP2-BLANK-SPIKE	48	1.00	2.00	2.00	100.0	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$