



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 7/6/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 14:14
In Date: 07/02/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:15
Out Date: 07/03/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB115194

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M2917-01	FK-BPM-04-017-ABCD	1	1.18	9.77	9.28	94.3	
M2917-01D	FK-BPM-04-017-ABCDDUP	2	1.17	9.77	9.3	94.5	
M2917-02	FK-BPM-04-017-ABCD	3	1.18	9.58	9.08	94.0	
M2918-01	FK-BPM-03-023-ABCD	4	1.16	9.61	8.84	90.9	
M2918-02	FK-BPM-03-023-ABCD	5	1.16	9.59	8.88	91.6	
M2919-01	FK-BPM-01-B9-B9A-B10-B10A	6	1.19	9.84	9.31	93.9	
M2919-02	FK-BPM-01-B9-B9A-B10-B10A	7	1.17	9.57	8.96	92.7	
M2922-01	SS-1 (0.5-1)	8	1.18	9.5	8.79	91.5	
M2922-02	SS-2 (0.5-1)	9	1.16	9.7	8.78	89.2	
M2922-03	SS-3 (0.5-1)	10	1.16	9.79	9.36	95.0	
M2922-04	SS-4 (0.5-1)	11	1.17	9.71	8.25	82.9	
M2922-05	SS-5 (0.5-1)	12	1.16	9.78	8.27	82.5	
M2924-05	SVOC-GPC-BLANK	13	1.00	2.00	2.00	100.0	blank
M2924-06	PEST-GPC-BLANK	14	1.00	2.00	2.00	100.0	blank
M2924-07	PEST-GPC-BLANK-SPIKE	15	1.00	2.00	2.00	100.0	blank
M2924-08	PCB-GPC-BLANK	16	1.00	2.00	2.00	100.0	blank
M2924-09	PCB-GPC-BLANK-SPIKE	17	1.00	2.00	2.00	100.0	blank
M2924-10	SVOC-GPC2-BLANK	18	1.00	2.00	2.00	100.0	blank
M2924-11	PEST-GPC2-BLANK	19	1.00	2.00	2.00	100.0	blank
M2924-12	PEST-GPC2-BLANK-SPIKE	20	1.00	2.00	2.00	100.0	blank
M2924-13	PCB-GCP2-BLANK	21	1.00	2.00	2.00	100.0	blank
M2924-14	PCB-GCP2-BLANK-SPIKE	22	1.00	2.00	2.00	100.0	blank
M2926-01	TP-1	28	1.17	9.93	9.1	90.5	
M2926-01D	TP-1DUP	29	1.17	9.78	9.12	92.3	
M2926-02	TP-1-EPH-2	30	1.16	9.85	8.92	89.3	
M2926-04	TP-2	31	1.16	9.69	8.78	89.3	
M2926-05	TP-2-EPH-3	32	1.18	9.41	8.57	89.8	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 7/6/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 14:14
In Date: 07/02/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:15
Out Date: 07/03/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB115194

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
M2926-07	TP-3	33	1.2	9.57	8.71	89.7	
M2926-08	TP-3-EPH-2	34	1.16	9.78	8.88	89.6	
M2931-01	HD-01-070221	37	1.18	9.86	8.86	88.5	
M2931-02	HD-01-070221-EPH-2	38	1.18	9.49	8.34	86.2	
M2932-01	DRUM-1	35	1.17	9.45	7.91	81.4	
M2932-02	DRUM-1	36	1.16	9.6	7.92	80.1	
M2933-01	1	23	1.17	9.54	8.57	88.4	
M2933-02	2	24	1.16	9.65	8.55	87.0	
M2933-03	3	25	1.17	9.29	8.24	87.1	
M2933-04	4	26	1.16	9.56	8.53	87.7	
M2933-05	5	27	1.16	9.67	8.54	86.7	
M2936-01	TR-05-070221	39	1.17	9.74	8.69	87.7	
M2936-02	TR-05-070221-EPH-2	40	1.18	9.8	8.88	89.3	
M2938-01	OK-01-070221	41	1.17	9.69	8.8	89.6	
M2938-02	OK-01-070221-EPH-2	42	1.19	9.77	8.83	89.0	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$