



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 7/9/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 07/08/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:16
Out Date: 07/09/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB115255

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M2969-01	SB01	1	1.15	9.76	8.55	85.9	
M2969-03	SB02	2	1.16	9.6	8.53	87.3	
M2969-05	SB03	3	1.11	9.7	8.27	83.4	
M2969-07	SB04	4	1.14	9.73	9.1	92.7	
M2969-09	SB05	5	1.15	9.55	8.6	88.7	
M2969-11	SB06	6	1.18	9.97	8.35	81.6	
M2969-13	SB07	7	1.15	9.88	8.82	87.9	
M2969-15	SB08	8	1.14	9.94	9.24	92.0	
M2969-15D	SB08DUP	9	1.16	9.96	9.26	92.0	
M2977-01	P001-B46C01-01	10	1.00	2.00	2.00	100.0	CONCRETE SAMPLE
M2981-01	T-1	11	1.14	9.74	8.76	88.6	
M2981-02	T-2	12	1.16	9.82	8.79	88.1	
M2981-03	T-3	13	1.11	9.85	8.77	87.6	
M2981-04	T-4	14	1.12	9.63	8.86	91.0	
M2981-05	T-5	15	1.13	9.8	8.71	87.4	
M2982-01	33301	16	1.00	2.00	2.00	100.0	WIPE SAMPLE
M2982-02	33303	17	1.00	2.00	2.00	100.0	WIPE SAMPLE
M2983-01	33298	18	1.00	2.00	2.00	100.0	WIPE SAMPLE
M2983-02	33299	19	1.00	2.00	2.00	100.0	WIPE SAMPLE
M2983-03	33300	20	1.00	2.00	2.00	100.0	WIPE SAMPLE
M2983-04	505	21	1.00	2.00	2.00	100.0	WIPE SAMPLE
M2983-05	506	22	1.00	2.00	2.00	100.0	WIPE SAMPLE
M2984-01	72-11940	23	1.15	9.66	9.05	92.8	
M2984-03	VNJ-228	24	1.16	9.85	9.61	97.2	
M2985-01	33302	25	1.00	2.00	2.00	100.0	WIPE SAMPLE
M2985-02	33304	26	1.00	2.00	2.00	100.0	WIPE SAMPLE
M2985-03	33305	27	1.00	2.00	2.00	100.0	WIPE SAMPLE
M2985-04	125	28	1.00	2.00	2.00	100.0	WIPE SAMPLE



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 7/9/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 07/08/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:16
Out Date: 07/09/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB115255

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
M2985-05	126	29	1.00	2.00	2.00	100.0	WIPE SAMPLE
M2985-06	127	30	1.00	2.00	2.00	100.0	WIPE SAMPLE
M2985-07	128	31	1.00	2.00	2.00	100.0	WIPE SAMPLE
M2987-01	CHRT-24368	32	1.14	9.95	2.11	11.0	SLUDGE SAMPLE
M2988-06	GRAB-1	33	1.16	9.78	8.7	87.5	
M2988-07	GRAB-2	34	1.11	9.53	8.72	90.4	
M2988-08	GRAB-3	35	1.13	9.61	8.69	89.2	
M2989-01	OR-03-070821	36	1.15	9.55	8.79	91.0	
M2989-02	OR-03-070821-EPH-2	37	1.16	9.58	9.11	94.4	
M2990-01	TP-D	38	1.14	9.94	8.64	85.2	
M2990-03	TP-D1	39	1.15	9.96	8.82	87.1	
M2990-04	TP-G	40	1.12	9.69	8.36	84.5	
M2990-06	TP-G1	41	1.16	9.96	8.79	86.7	
M2990-07	TP-K	42	1.15	9.61	9.22	95.4	
M2990-09	TP-K1	43	1.16	9.6	9.24	95.7	
M2990-10	TP-I	44	1.13	9.85	9.06	90.9	
M2990-12	TP-I1	45	1.1	9.67	8.94	91.5	
M2990-13	TP-J	46	1.19	9.68	9.33	95.9	
M2990-15	TP-J1	47	1.16	9.94	9.58	95.9	
M2990-16	TP-E	48	1.11	9.84	9.23	93.0	
M2990-18	TP-E1	49	1.13	9.75	9.08	92.2	
M2991-01	33310	50	1.00	2.00	2.00	100.0	WIPE SAMPLE
M2991-02	33311	51	1.00	2.00	2.00	100.0	WIPE SAMPLE
M2991-03	33312	52	1.00	2.00	2.00	100.0	WIPE SAMPLE

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$