



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 7/15/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 17:10
In Date: 07/14/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:15
Out Date: 07/15/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB115355

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M3031-01	HR-03-071421	7	1.14	9.77	8.89	89.8	
M3031-02	HR-03-071421-EPH-2	8	1.16	9.91	9.04	90.1	
M3032-01	NB-07-071421	9	1.18	9.68	9.29	95.4	
M3032-02	NB-07-071421-EPH-2	10	1.16	9.84	9.44	95.4	
M3037-01	DS-1	1	1.13	9.52	8.73	90.6	
M3037-02	DN-1	2	1.19	9.63	8.97	92.2	
M3037-03	CS-1	3	1.14	9.94	8.86	87.7	
M3037-04	CN-1	4	1.18	9.73	8.48	85.4	
M3037-05	BN-1	5	1.19	9.74	8.85	89.6	
M3037-06	AN-1	6	1.19	9.9	8.74	86.7	
M3038-01	SB-06-5.5-6.0	11	1.17	9.7	8.94	91.1	
M3038-02	SB-06-COMP	12	1.15	9.95	9.33	93.0	
M3038-03	SB-07-5.5-6.0	13	1.15	9.58	8.57	88.0	
M3038-04	SB-07-COMP	14	1.14	9.83	8.82	88.4	
M3038-05	SB-09-5.5-6.0	15	1.15	9.63	8.1	82.0	
M3038-06	SB-09-COMP	16	1.12	9.68	8.44	85.5	
M3038-07	SB-11-5.5-6.0	17	1.15	9.95	8.55	84.1	
M3038-08	SB-11-COMP	18	1.11	9.76	8.34	83.6	
M3038-09	SB-10-5.5-6.0	19	1.15	9.96	9.5	94.8	
M3038-10	SB-10-COMP	20	1.14	9.85	9.41	94.9	
M3039-01	SB-12-5.5-6.0	21	1.15	9.95	9.26	92.2	
M3039-02	SB-12-COMP	22	1.12	9.59	8.76	90.2	
M3039-03	SB-01-9.5-10.0	23	1.14	9.98	8.7	85.5	
M3039-04	SB-01-COMP	24	1.19	9.82	8.07	79.7	
M3039-05	SB-02-9.0-9.5	25	1.19	9.81	8.07	79.8	
M3039-06	SB-02-COMP	26	1.1	9.62	8.44	86.2	
M3039-07	SB-03-9.5-10.0	27	1.14	9.94	9.77	98.1	
M3039-08	SB-03-COMP	28	1.18	9.81	9.56	97.1	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 7/15/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 17:10
In Date: 07/14/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:15
Out Date: 07/15/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB115355

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
M3039-09	SB-04-8.0-8.5	29	1.13	9.94	9.14	90.9	
M3039-10	SB-04-COMP	30	1.16	9.58	8.82	91.0	
M3039-11	SB-05-11.5-12.0	31	1.15	9.7	8.7	88.3	
M3039-12	SB-05-COMP	32	1.16	9.98	9.12	90.2	
M3039-13	SB-08-11.5-12.0	33	1.1	9.71	9.06	92.5	
M3039-14	SB-08-COMP	34	1.19	9.79	8.83	88.8	
M3040-01	FK-PS-01-076	35	1.15	9.5	7.9	80.8	
M3040-02	FK-PS-01-076	36	1.14	9.58	7.95	80.7	
M3043-01	FK-PS-01-077	37	1.16	9.62	8.01	81.0	
M3043-02	FK-PS-01-077	38	1.16	9.69	8.3	83.7	
M3046-01	FK-EF-01-019-ABCD	39	1.13	9.72	9.04	92.1	
M3046-02	FK-EF-01-019-ABCD	40	1.15	9.54	8.89	92.3	
M3047-01	FK-EF-020-ABCD	41	1.16	9.73	9.1	92.6	
M3047-02	FK-EF-020-ABCD	42	1.16	9.86	9.16	92.0	
M3048-01	FK-BPM-01-AA1-AA1A-BB1-BB1A	43	1.12	9.74	9.13	92.9	
M3048-02	FK-BPM-01-AA1-AA1A-BB1-BB1A	44	1.16	9.74	9.22	93.9	
M3049-01	FK-BPM-03-024-ABCD	45	1.16	9.79	9.1	92.0	
M3049-02	FK-BPM-03-024-ABCD	46	1.18	9.72	9.27	94.7	
M3050-01	FK-BPM-04-018-ABCD	47	1.13	9.51	9.05	94.5	
M3050-01D	FK-BPM-04-018-ABCD DUP	48	1.16	9.53	9.06	94.4	
M3050-02	FK-BPM-04-018-ABCD	49	1.17	9.97	9.28	92.2	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$