



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 7/16/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 17:10
In Date: 07/15/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:10
Out Date: 07/16/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB115374

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M3051-01	41260	30	1.18	9.74	8.62	86.9	
M3051-02	30778	31	1.00	2.00	2.00	100.0	OIL SAMPLE
M3051-04	41904	32	1.00	2.00	2.00	100.0	OIL SAMPLE
M3051-06	41272	33	1.00	2.00	2.00	100.0	OIL SAMPLE
M3051-08	41213	34	1.14	9.98	9.00	88.9	sludge
M3052-01	BUR-1085	1	1.00	2.00	2.00	100.0	WIPE SAMPLE
M3052-02	BUR-1086	2	1.00	2.00	2.00	100.0	WIPE SAMPLE
M3052-03	BUR-1087	3	1.00	2.00	2.00	100.0	WIPE SAMPLE
M3053-01	BC204799LLK-END-1-1	22	1.00	2.00	2.00	100.0	PILC REEL
M3053-02	BC204799LLK-END-1-2	23	1.00	2.00	2.00	100.0	PILC REEL
M3053-03	BC204799LLK-END-2-1	24	1.00	2.00	2.00	100.0	PILC REEL
M3053-04	BC204799LLK-END-2-2	25	1.00	2.00	2.00	100.0	PILC REEL
M3053-05	BL226480L0Z-END-1-1	26	1.00	2.00	2.00	100.0	PILC REEL
M3053-06	BL226480L0Z-END-1-2	27	1.00	2.00	2.00	100.0	PILC REEL
M3053-07	HEH768H-END-1-1	28	1.00	2.00	2.00	100.0	PILC REEL
M3053-08	HEH768H-END-1-2	29	1.00	2.00	2.00	100.0	PILC REEL
M3054-01	HIB135S-END-1-1	4	1.00	2.00	2.00	100.0	PILC REEL
M3054-02	HIB135S-END-1-2	5	1.00	2.00	2.00	100.0	PILC REEL
M3054-03	FCA568Y-END-1-1	6	1.00	2.00	2.00	100.0	PILC REEL
M3054-04	FCA568Y-END-1-2	7	1.00	2.00	2.00	100.0	PILC REEL
M3054-05	FCA568Y-END-1-3	8	1.00	2.00	2.00	100.0	PILC REEL
M3054-06	FCA568Y-END-1-4	9	1.00	2.00	2.00	100.0	PILC REEL
M3054-07	HIA869E-END-1-1	10	1.00	2.00	2.00	100.0	PILC REEL
M3054-08	HIA869E-END-1-2	11	1.00	2.00	2.00	100.0	PILC REEL
M3054-09	HIA869E-END-1-3	12	1.00	2.00	2.00	100.0	PILC REEL
M3054-10	HIA869E-END-1-4	13	1.00	2.00	2.00	100.0	PILC REEL
M3055-01	KMA6125-END1-1	14	1.00	2.00	2.00	100.0	PILC REEL
M3055-02	KMA6125-END1-1	15	1.00	2.00	2.00	100.0	PILC REEL



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 7/16/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 17:10
In Date: 07/15/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:10
Out Date: 07/16/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB115374

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M3056-01	RO-VNJ-239	16	1.18	9.68	8.61	87.4	
M3056-02	RO-VNJ-239-EPH-2	17	1.16	9.78	8.77	88.3	
M3061-01	PB-2021-WW-000-08	18	1.15	9.97	6.35	59.0	
M3065-01	TR-04-071521	19	1.15	9.58	9.06	93.8	
M3065-01D	TR-04-071521DUP	20	1.17	9.6	9.08	93.8	
M3065-02	TR-04-071521-EPH-2	21	1.16	9.71	9.18	93.8	
M3067-01	RW-3	35	1.00	2.00	2.00	100.0	OIL SAMPLE

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$