



PERCENT SOLID

Supervisor: apatel  
Analyst: jignesh  
Date: 7/21/2021

OVENTEMP IN Celsius(°C): 107  
Time IN: 17:15  
In Date: 07/20/2021  
Weight Check 1.0g: 1.00  
Weight Check 10g: 10.00  
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103  
Time OUT: 08:28  
Out Date: 07/21/2021  
Weight Check 1.0g: 1.00  
Weight Check 10g: 10.00  
BalanceID: M SC-4  
Thermometer ID: %SOLIDS-OVEN

QC:LB115442

| Lab ID    | Client SampleID        | Dish # | Dish Wt(g) (A) | Dish + Sample Wt(g) (B) | Dish+Dry Sample Wt(g) (C) | % Solid | Comments        |
|-----------|------------------------|--------|----------------|-------------------------|---------------------------|---------|-----------------|
| M3059-03  | 01-A-01-B-01-C         | 15     | 1.00           | 2.00                    | 2.00                      | 100.0   | CAULKING sample |
| M3059-04  | 02-A-02-B-02-C         | 16     | 1.00           | 2.00                    | 2.00                      | 100.0   | CAULKING sample |
| M3098-01  | F0604S                 | 1      | 1.12           | 9.87                    | 4.95                      | 43.8    | SLUDGE SAMPLE   |
| M3099-01  | BC245062L0Z-END-1-1    | 2      | 1.00           | 2.00                    | 2.00                      | 100.0   | PILC REEL       |
| M3099-02  | BC245062L0Z-END-1-2    | 3      | 1.00           | 2.00                    | 2.00                      | 100.0   | PILC REEL       |
| M3099-03  | BC123196LOK-END-1-1    | 4      | 1.00           | 2.00                    | 2.00                      | 100.0   | PILC REEL       |
| M3099-04  | BC123196LOK-END-1-2    | 5      | 1.00           | 2.00                    | 2.00                      | 100.0   | PILC REEL       |
| M3100-01  | BC260356LOK-END-1-1    | 6      | 1.00           | 2.00                    | 2.00                      | 100.0   | PILC REEL       |
| M3100-02  | BC260356LOK-END-1-2    | 7      | 1.00           | 2.00                    | 2.00                      | 100.0   | PILC REEL       |
| M3100-03  | BC260053LOK-END-1-1    | 8      | 1.00           | 2.00                    | 2.00                      | 100.0   | PILC REEL       |
| M3100-04  | BC260053LOK-END-1-2    | 9      | 1.00           | 2.00                    | 2.00                      | 100.0   | PILC REEL       |
| M3100-05  | JEB471U-END-1-1        | 10     | 1.00           | 2.00                    | 2.00                      | 100.0   | PILC REEL       |
| M3100-06  | JEB471U-END-1-2        | 11     | 1.00           | 2.00                    | 2.00                      | 100.0   | PILC REEL       |
| M3101-01  | 33324                  | 12     | 1.00           | 2.00                    | 2.00                      | 100.0   | OILY DEBRIS     |
| M3109-01  | PX-02D                 | 13     | 1.15           | 9.96                    | 8.48                      | 83.2    |                 |
| M3109-01D | PX-02DDUP              | 14     | 1.14           | 9.98                    | 8.5                       | 83.3    |                 |
| M3114-01  | 20210541-COM-14        | 17     | 1.16           | 9.96                    | 4.88                      | 42.3    | sediment sample |
| M3114-02  | 20210550-COM-23        | 18     | 1.17           | 9.77                    | 4.54                      | 39.2    | sediment sample |
| M3114-03  | 20210559-COM-32        | 19     | 1.17           | 9.69                    | 4.66                      | 41.0    | sediment sample |
| M3114-04  | 20210585-COM-8         | 20     | 1.13           | 9.98                    | 5.38                      | 48.0    | sediment sample |
| M3114-05  | 20210594-AMENDED-COM-8 | 21     | 1.16           | 9.89                    | 5.42                      | 48.8    | sediment sample |

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$