



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 8/13/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 14:05
In Date: 07/27/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:13
Out Date: 07/28/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB115551

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
M3187-01	GB5Q7	1	1.11	9.75	9.55	97.7	
M3187-02	GB5Q8	2	1.16	9.59	9.4	97.7	
M3187-03	GB5Q9	3	1.16	9.81	9.73	99.1	
M3187-04	GB5R0	4	1.12	9.82	9.75	99.2	
M3187-05	GB5R1	5	1.13	9.73	8.86	89.9	
M3187-06	GB5R1D	6	1.13	9.73	8.86	89.9	
M3187-07	GB5R1S	7	1.13	9.73	8.86	89.9	
M3187-08	GB5R2	8	1.15	9.52	9.25	96.8	
M3187-09	GB5R3	9	1.18	9.65	9.48	98.0	
M3187-10	GB5R4	10	1.15	9.59	9.37	97.4	
M3187-11	GB5R5	11	1.13	9.9	9.62	96.8	
M3187-12	GB5R6	12	1.16	9.74	9.62	98.6	
M3187-13	GB5R7	13	1.11	9.86	9.75	98.7	
M3187-14	GB5R8	14	1.14	9.94	9.5	95.0	
M3187-15	GB5R9	15	1.13	9.71	9.65	99.3	
M3187-16	GB5S0	16	1.16	9.88	9.64	97.2	
M3187-17	GB5S1	17	1.15	9.9	9.53	95.8	
M3187-18	GB5S2	18	1.18	9.51	8.99	93.8	
M3187-19	GB5S3	19	1.00	2.00	2.00	100.0	P.E. SAMPLE

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$