



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 8/24/2021

OVENTEMP IN Celsius(°C): 108
Time IN: 16:45
In Date: 08/11/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:00
Out Date: 08/12/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB115789

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
M3359-01	JEH585T-END-1-1	34	1.00	2.00	2.00	100.0	PILC REEL
M3359-02	JEH585T-END-1-2	35	1.00	2.00	2.00	100.0	PILC REEL
M3359-03	NO-NUMBER-END-1-1	36	1.00	2.00	2.00	100.0	PILC REEL
M3359-04	NO-NUMBER-END-1-2	37	1.00	2.00	2.00	100.0	PILC REEL
M3359-05	FDL007D-END-1-1	38	1.00	2.00	2.00	100.0	PILC REEL
M3359-06	FDL007D-END-1-2	39	1.00	2.00	2.00	100.0	PILC REEL
M3359-07	HED160B-END-1-1	40	1.00	2.00	2.00	100.0	PILC REEL
M3359-08	HED160B-END-1-2	41	1.00	2.00	2.00	100.0	PILC REEL
M3361-01	41944	51	1.13	9.77	9.7	99.2	
M3361-02	41945	52	1.14	9.8	9.76	99.5	
M3361-03	41946	53	1.19	9.69	9.67	99.8	
M3361-04	41947	54	1.12	9.63	9.61	99.8	
M3361-05	41910	55	1.00	2.00	2.00	100.0	OIL SAMPLE
M3361-06	44-45-46-47	56	1.13	9.68	9.65	99.6	
M3363-01	PCB-080521-TB-01	57	1.00	2.00	2.00	100.0	CAULKING SAMPLE
M3363-02	PCB-080521-TB-02	58	1.00	2.00	2.00	100.0	CAULKING SAMPLE
M3363-03	PCB-080521-TB-03	59	1.00	2.00	2.00	100.0	CAULKING SAMPLE
M3368-01	SWCS-01-0002	1	1.15	9.64	9.16	94.3	
M3368-02	SWCS-01-0205	2	1.18	9.53	8.99	93.5	
M3368-03	SWCS-01-0510	3	1.16	9.94	8.84	87.5	
M3368-04	SWCS-01-1015	4	1.14	9.94	8.12	79.3	
M3368-05	SWCS-02-0002	5	1.14	9.94	9.6	96.1	
M3368-06	SWCS-02-0205	6	1.15	9.86	8.78	87.6	
M3368-07	SWCS-02-0510	7	1.16	9.94	7.75	75.1	
M3368-08	SWCS-07-0002	8	1.17	9.7	8.76	89.0	
M3368-09	SWCS-07-0205	9	1.12	9.68	8.55	86.8	
M3368-10	SWCS-07-0510	10	1.16	9.85	8.62	85.8	
M3368-11	SWCS-07-1015	11	1.16	9.9	9.14	91.3	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 8/24/2021

OVENTEMP IN Celsius(°C): 108
Time IN: 16:45
In Date: 08/11/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:00
Out Date: 08/12/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB115789

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M3368-12	SWCS-06-0002	12	1.15	9.61	8.89	91.5	
M3368-13	SWCS-06-0205	13	1.13	9.71	8.59	86.9	
M3368-14	SWCS-06-0510	14	1.15	9.75	8.36	83.8	
M3368-15	SWCS-06-1015	15	1.13	9.73	8.7	88.0	
M3368-16	SWCS-05-0002	16	1.16	9.54	8.9	92.4	
M3368-17	2021-0810-FD-01	17	1.19	9.7	8.36	84.3	
M3368-18	SWCS-05-0205	18	1.1	9.86	9.19	92.4	
M3368-19	SWCS-05-0510	19	1.19	9.82	8.6	85.9	
M3368-20	SWCS-04-0002	20	1.16	9.5	9.11	95.3	
M3368-21	SWCS-04-0205	21	1.1	9.91	8.84	87.9	
M3368-22	SWCS-04-0510	22	1.16	9.81	7.82	77.0	
M3371-01	JET-GROUT	50	1.18	9.97	6.5	60.5	
M3372-01	HKA318J-END-1-1	42	1.00	2.00	2.00	100.0	PILC REEL
M3372-02	HKA318J-END-1-2	43	1.00	2.00	2.00	100.0	PILC REEL
M3372-03	KEA083T-END-1-1	44	1.00	2.00	2.00	100.0	PILC REEL
M3372-04	KEA083T-END-1-2	45	1.00	2.00	2.00	100.0	PILC REEL
M3372-05	KMH498D-END-1-1	46	1.00	2.00	2.00	100.0	PILC REEL
M3372-06	KMH498D-END-1-2	47	1.00	2.00	2.00	100.0	PILC REEL
M3372-07	KMH498D-END-2-1	48	1.00	2.00	2.00	100.0	PILC REEL
M3372-08	KMH498D-END-2-2	49	1.00	2.00	2.00	100.0	PILC REEL
M3375-01	T-1	23	1.14	9.86	8.65	86.1	
M3375-01D	T-1DUP	24	1.16	9.89	8.67	86.0	
M3375-02	T-2	25	1.12	9.95	8.89	88.0	
M3375-03	T-3	26	1.16	9.72	8.71	88.2	
M3375-04	T-4	27	1.14	9.77	8.9	89.9	
M3375-05	T-5	28	1.13	9.53	8.46	87.3	
M3375-06	T-6	29	1.16	9.78	9.09	92.0	
M3375-07	T-7	30	1.16	9.91	8.95	89.0	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 8/24/2021

OVENTEMP IN Celsius(°C): 108
Time IN: 16:45
In Date: 08/11/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:00
Out Date: 08/12/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB115789

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M3375-08	T-8	31	1.17	9.6	8.82	90.7	
M3375-09	P-1	32	1.15	9.97	9.06	89.7	
M3375-10	P-2	33	1.16	9.98	9.47	94.2	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$