



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 8/13/2021

OVENTEMP IN Celsius(°C): 108
Time IN: 17:40
In Date: 08/12/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 104
Time OUT: 08:25
Out Date: 08/13/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB115812

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M3369-01	S2-04-B	1	1.15	9.63	8.63	88.2	
M3384-01	MAIN-BODY	2	1.00	2.00	2.00	100.0	OIL SAMPLE
M3384-02	TAP-CHANGER	3	1.00	2.00	2.00	100.0	OIL SAMPLE
M3385-01	MAIN-BODY	4	1.00	2.00	2.00	100.0	OIL SAMPLE
M3385-02	TAP-CHANGER	5	1.00	2.00	2.00	100.0	OIL SAMPLE
M3385-03	CAPACIFORMER	6	1.00	2.00	2.00	100.0	OIL SAMPLE
M3385-04	PRIMARY	7	1.00	2.00	2.00	100.0	OIL SAMPLE
M3386-01	TRANSFO-MAIN	8	1.00	2.00	2.00	100.0	OIL SAMPLE
M3386-02	TRANSFO-CAP	9	1.00	2.00	2.00	100.0	OIL SAMPLE
M3386-03	TAP	10	1.00	2.00	2.00	100.0	OIL SAMPLE
M3386-04	PRIMARY	11	1.00	2.00	2.00	100.0	OIL SAMPLE
M3387-01	MAIN	12	1.00	2.00	2.00	100.0	OIL SAMPLE
M3387-02	POT	13	1.00	2.00	2.00	100.0	OIL SAMPLE
M3387-03	TAP	14	1.00	2.00	2.00	100.0	OIL SAMPLE
M3387-04	CAP	15	1.00	2.00	2.00	100.0	OIL SAMPLE
M3388-01	VNJ-215	16	1.16	9.58	8.22	83.8	
M3388-03	VNJ-215-EPH-2	17	1.16	9.58	8.22	83.8	
M3389-01	CRT-21	18	1.19	9.73	8.99	91.3	
M3389-03	CRT-21-EPH-2	19	1.18	9.62	9.04	93.1	
M3389-06	CRT-18	20	1.19	9.72	9.14	93.2	
M3389-08	CRT-18-EPH-2	21	1.16	9.64	9.06	93.2	
M3389-11	CRT-17	22	1.11	9.98	9.27	92.0	
M3389-13	CRT-17-EPH-2	23	1.17	9.72	8.86	89.9	
M3390-01	SWCS-03-0002	42	1.17	9.97	9.4	93.5	
M3390-02	SWCS-03-0205	43	1.17	9.62	8.87	91.1	
M3390-03	SWCS-03-0510	44	1.17	9.81	7.12	68.9	
M3390-04	SWCS-03-1015	45	1.17	9.98	9.28	92.1	
M3390-05	M3390-04MS	46	1.17	9.98	9.28	92.1	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 8/13/2021

OVENTEMP IN Celsius(°C): 108
Time IN: 17:40
In Date: 08/12/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 104
Time OUT: 08:25
Out Date: 08/13/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB115812

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
M3390-06	M3390-04MSD	47	1.17	9.98	9.28	92.1	
M3390-07	SWCS-08-0002	48	1.17	9.98	9.31	92.4	
M3390-08	SWCS-08-0205	49	1.14	9.84	9.18	92.4	
M3390-09	SWCS-08-0510	50	1.13	9.96	9.24	91.8	
M3390-10	SWCS-08--1015	51	1.16	9.93	9.17	91.3	
M3390-11	SWCS-09-0002	52	1.13	9.82	8.15	80.8	
M3390-12	SWCS-09-0205	53	1.16	9.66	8.16	82.4	
M3390-13	SWCS-09-0510	54	1.14	9.92	9.3	92.9	
M3390-14	2021-08-11-FD-02	55	1.15	9.98	9.6	95.7	
M3390-15	SWCS-09-1015	56	1.16	9.89	8.54	84.5	
M3390-16	SWCS-10-0002	57	1.14	9.56	8.67	89.4	
M3390-17	SWCS-10-0205	58	1.1	9.7	8.39	84.8	
M3390-18	SWCS-10-0510	59	1.15	9.97	9.22	91.5	
M3390-19	M3390-18MS	60	1.15	9.97	9.22	91.5	
M3390-20	M3390-18MSD	61	1.15	9.97	9.22	91.5	
M3390-21	SWCS-10-1015	62	1.19	9.65	8.44	85.7	
M3396-01	PIPE-A	24	1.00	2.00	2.00	100.0	WIPE SAMPLE
M3396-02	PIPE-B	25	1.00	2.00	2.00	100.0	WIPE SAMPLE
M3396-03	PIPE-C	26	1.00	2.00	2.00	100.0	WIPE SAMPLE
M3396-04	PIPE-D	27	1.00	2.00	2.00	100.0	WIPE SAMPLE
M3396-05	PIPE-E	28	1.00	2.00	2.00	100.0	WIPE SAMPLE
M3396-06	PIPE-F	29	1.00	2.00	2.00	100.0	WIPE SAMPLE
M3401-01	BR-2-081221	30	1.16	9.77	9.1	92.2	
M3403-01	VNJ-213	31	1.00	2.00	2.00	100.0	CONCRETE SAMPLE
M3404-01	TP1-1-6FT	32	1.14	9.94	8.9	88.2	
M3404-02	TP-1-6FT-EPH-2	33	1.14	9.94	8.9	88.2	
M3404-04	TP2-1-6FT	34	1.17	9.86	8.85	88.4	
M3404-05	TP2-1-6FT-EPH-2	35	1.17	9.86	8.85	88.4	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 8/13/2021

OVENTEMP IN Celsius(°C): 108
Time IN: 17:40
In Date: 08/12/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 104
Time OUT: 08:25
Out Date: 08/13/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB115812

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M3404-07	TP3-1-6FT	36	1.19	9.98	8.95	88.3	
M3404-08	TP3-1-6FT-EPH-2	37	1.19	9.98	8.95	88.3	
M3404-10	CONCRETE-SLAB	38	1.00	2.00	2.00	100.0	CONCRETE SAMPLE
M3405-01	CL-01-081221	39	1.19	9.8	9.00	90.7	
M3405-02	CL-01-081221-EPH-2	40	1.19	9.92	9.22	92.0	
M3405-02D	CL-01-081221-EPH-2DUP	41	1.17	9.94	9.24	92.0	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$