



PERCENT SOLID

Supervisor: ketankumar
Analyst: jignesh
Date: 8/19/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 14:00
In Date: 08/18/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:11
Out Date: 08/19/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB115912

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
M3365-01	DBK80	1	1.15	9.9	8.05	78.9	
M3365-02	DBK81	2	1.16	9.75	7.59	74.9	
M3365-03	DBK83	3	1.13	9.85	7.98	78.6	
M3365-04	DBK85	4	1.15	9.95	7.78	75.3	
M3365-05	DBK87	5	1.13	9.83	7.19	69.7	
M3365-06	DBK89	6	1.11	9.71	7.17	70.5	
M3365-07	DBK91	7	1.16	9.96	8.6	84.5	
M3365-08	DBKA0	8	1.16	9.98	8.52	83.4	
M3365-09	VHBLK001	9	1.00	2.00	2.00	100.0	VHBLK
M3365-10	DBKB4	10	1.14	9.96	7.37	70.6	
M3365-11	DBKB5	11	1.1	9.68	6.91	67.7	
M3365-12	DBKB6	12	1.18	9.69	6.77	65.7	
M3365-13	DBKB8	13	1.16	9.97	7.21	68.7	
M3365-14	DBKC0	14	1.16	9.97	8.12	79.0	
M3365-15	DBKCOMS	15	1.16	9.97	8.12	79.0	
M3365-16	DBKCOMSD	16	1.16	9.97	8.12	79.0	
M3365-17	DBK76	17	1.00	2.00	2.00	100.0	P.E. SAMPLE
M3365-18	DBK77	18	1.00	2.00	2.00	100.0	P.E. SAMPLE

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$