



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 8/20/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 15:35
In Date: 08/19/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-01

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:11
Out Date: 08/20/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB115935

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
M3472-01	C0AB9	1	1.12	9.55	7.78	79.0	
M3472-02	C0AC0	2	1.12	9.84	8.74	87.4	
M3472-03	C0AC1	3	1.16	9.67	7.61	75.8	
M3472-04	C0AC2	4	1.18	9.72	8.64	87.4	
M3472-05	C0AS3	5	1.16	9.81	7.69	75.5	
M3472-06	C0AX0	6	1.17	9.74	8.13	81.2	
M3472-07	C0B89	7	1.12	9.55	8.81	91.2	
M3472-09	C0AB9	8	1.12	9.55	7.78	79.0	
M3472-10	C0AC0	9	1.12	9.84	8.74	87.4	
M3472-11	C0AC1	10	1.16	9.67	7.61	75.8	
M3472-12	C0AC2	11	1.18	9.72	8.64	87.4	
M3472-13	C0AS3	12	1.16	9.81	7.69	75.5	
M3472-14	C0AX0	13	1.17	9.74	8.13	81.2	
M3472-19	C0B89	14	1.14	9.75	8.81	89.1	
M3472-20	VHBLK001	15	1.00	2.00	2.00	100.0	VHBLK

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$