



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 8/23/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 15:00
In Date: 08/21/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-01

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:22
Out Date: 08/22/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB115984

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
M3486-01	C0AR6	1	1.12	9.52	7.8	79.5	
M3486-02	C0AW2	2	1.16	9.74	8.84	89.5	
M3486-03	C0AW9	3	1.18	9.64	8.12	82.0	
M3486-04	C0AB7	4	1.14	9.54	8.03	82.0	
M3486-05	C0AB8	5	1.15	9.97	8.7	85.6	
M3486-06	C0AR5	6	1.18	9.66	7.53	74.9	
M3486-07	C0AR9	7	1.12	9.86	7.51	73.1	
M3486-08	C0AW0	8	1.13	9.55	8.11	82.9	
M3486-09	C0AW5	9	1.15	9.96	8.08	78.7	
M3486-10	VHBLK001	10	1.00	2.00	2.00	100.0	VHBLK
M3486-11	C0AY8	11	1.15	9.83	8.01	79.0	
M3486-12	C0AY9	12	1.16	9.94	8.34	81.8	
M3486-13	C0AZ4	13	1.15	9.95	6.87	65.0	
M3486-14	C0B05	14	1.1	9.86	8.15	80.5	
M3486-15	C0B19	15	1.12	9.76	8.08	80.6	
M3486-16	C0B24	16	1.13	9.63	7.55	75.5	
M3486-17	C0B25	17	1.12	9.98	8.19	79.8	
M3486-18	C0B30	18	1.18	9.67	8.77	89.4	
M3486-19	C0B41	19	1.18	9.6	8.3	84.6	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$