



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 8/27/2021

OVENTEMP IN Celsius(°C): 108
Time IN: 17:00
In Date: 08/26/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 104
Time OUT: 08:16
Out Date: 08/27/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB116030

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M3541-01	1	1	1.14	9.97	8.25	80.5	
M3541-02	2	2	1.14	9.94	8.3	81.4	
M3541-03	3	3	1.12	9.92	7.89	76.9	
M3541-04	4	4	1.17	9.89	7.94	77.6	
M3541-05	5	5	1.16	9.71	7.87	78.5	
M3545-01	NCS-COMP-1	6	1.16	9.63	8.97	92.2	
M3545-02	NCS-EPH-2	7	1.16	9.64	8.92	91.5	
M3545-03	NCS-COMP-2	8	1.17	9.74	9.04	91.8	
M3545-04	NCS-EPH-4	9	1.15	9.96	8.71	85.8	
M3546-01	TPR-4-COMP	10	1.15	9.75	8.22	82.2	
M3546-02	TRP-4-EPH-2	11	1.14	9.91	8.77	87.0	
M3546-03	TRP-7-COMP	12	1.16	9.97	8.74	86.0	
M3546-04	TRP-7-EPH-2	13	1.11	9.98	8.44	82.6	
M3547-01	PE-1	14	1.19	9.93	9.04	89.8	
M3547-02	PE-1D	15	1.17	9.57	8.92	92.3	
M3547-03	PE-2	16	1.16	9.98	8.79	86.5	
M3547-04	PE-3	17	1.19	9.71	8.7	88.1	
M3547-05	PE-4	18	1.1	9.88	9.17	91.9	
M3547-06	PE-5	19	1.18	9.77	8.46	84.7	
M3547-07	PE-6	20	1.16	9.83	8.56	85.4	
M3547-08	PE-7	21	1.17	9.66	9.00	92.2	
M3547-08D	PE-7DUP	22	1.19	9.69	9.01	92.0	
M3548-01	VNJ-235	23	1.14	9.77	9.02	91.3	
M3548-02	VNJ-235-EPH-2	24	1.14	9.56	8.83	91.3	
M3550-01	WC-1	33	1.13	9.84	8.55	85.2	
M3551-01	250EM-139-0002	27	1.13	9.62	8.8	90.3	
M3551-02	250EM-139-0204	28	1.16	9.68	7.73	77.1	
M3551-03	250EM-140-0002	29	1.17	9.96	9.33	92.8	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 8/27/2021

OVENTEMP IN Celsius(°C): 108
Time IN: 17:00
In Date: 08/26/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 104
Time OUT: 08:16
Out Date: 08/27/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB116030

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M3551-04	250EM-140-0204	30	1.16	9.9	9.45	94.9	
M3557-01	OK-02-082621	25	1.12	9.59	9.03	93.4	
M3557-02	OK-02-082621-EPH-2	26	1.16	9.65	9.11	93.6	
M3558-01	TPR-14-COMP	31	1.17	9.88	8.7	86.5	
M3558-02	TPR-14-EPH-2	32	1.15	9.89	8.89	88.6	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$