



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 9/10/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 17:20
In Date: 09/09/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:17
Out Date: 09/10/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB116213

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M3682-01	WC-2	1	1.14	9.78	8.78	88.4	
M3682-02	WC-3	2	1.16	9.67	8.26	83.4	
M3682-03	WC-4	3	1.16	9.96	8.74	86.1	
M3683-01	01-A-01-B-01-C	6	1.00	2.00	2.00	100.0	CAULKING SAMPLE
M3684-01	A018015	43	1.16	9.97	8.81	86.8	
M3684-02	A018145	44	1.15	9.59	9.00	93.0	
M3684-03	DUP-210908A	45	1.14	9.77	9.18	93.2	
M3684-04	F003130	46	1.15	9.69	8.55	86.7	
M3684-05	DUP-210908B	47	1.14	9.77	8.61	86.6	
M3684-06	F004130	48	1.15	9.97	8.91	88.0	
M3684-07	A019015	49	1.14	9.53	8.76	90.8	
M3684-08	A019145	50	1.14	9.96	9.6	95.9	
M3684-09	F002050	51	1.14	9.56	8.97	93.0	
M3684-10	A017015	52	1.14	9.83	8.3	82.4	
M3684-11	A017145	53	1.14	9.77	8.52	85.5	
M3684-12	F005130	54	1.14	9.72	8.45	85.2	
M3684-13	F006130	55	1.13	9.67	8.85	90.4	
M3684-14	A008015	56	1.14	9.66	8.82	90.1	
M3684-15	A008145	57	1.15	9.97	9.57	95.5	
M3684-16	A007015	58	1.16	9.96	9.24	91.8	
M3684-17	A007145	59	1.14	9.74	9.13	92.9	
M3684-18	A001015	60	1.16	9.95	8.92	88.3	
M3684-19	A001145	61	1.15	9.51	8.22	84.6	
M3684-20	A006015	62	1.14	9.95	9.63	96.4	
M3684-21	A006070	63	1.14	9.97	9.42	93.8	
M3685-01	327	38	1.14	9.96	8.99	89.0	
M3685-02	327-EPH-2	39	1.14	9.7	9.12	93.2	
M3686-01	MAIN-GUARD-HOUSE-JET-GR OUT	4	1.15	9.57	6.5	63.5	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 9/10/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 17:20
In Date: 09/09/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:17
Out Date: 09/10/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB116213

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M3686-01D	MAIN-GUARD-HOUSE-JET-GR OUTDUP	5	1.17	9.59	6.53	63.7	
M3689-01	256	40	1.00	2.00	2.00	100.0	OILY-DEBRIS
M3691-01	OK-03-090921	41	1.12	9.8	8.91	89.7	
M3691-02	OK-03-090921-EPH-2	42	1.19	9.56	8.73	90.1	
M3692-01	WM-1	7	1.15	9.95	8.94	88.5	
M3692-02	WM-1A	8	1.16	9.98	9.19	91.0	
M3693-01	WFM-1	9	1.17	9.86	8.96	89.6	
M3694-01	S-1	10	1.15	9.96	8.49	83.3	
M3694-02	S-2	11	1.14	9.94	8.86	87.7	
M3694-03	S-3	12	1.1	9.87	8.66	86.2	
M3694-04	S-4	13	1.12	9.95	9.32	92.9	
M3696-01	HIH012E-END-1-1	14	1.00	2.00	2.00	100.0	PILC REEL
M3696-02	HIH012E-END-1-2	15	1.00	2.00	2.00	100.0	PILC REEL
M3696-03	HIH012E-END-2-1	16	1.00	2.00	2.00	100.0	PILC REEL
M3696-04	HIH012E-END-2-2	17	1.00	2.00	2.00	100.0	PILC REEL
M3696-05	HIH012E-END-3-1	18	1.00	2.00	2.00	100.0	PILC REEL
M3696-06	HIH012E-END-3-2	19	1.00	2.00	2.00	100.0	PILC REEL
M3696-07	869141B-3B-END-1-1	20	1.00	2.00	2.00	100.0	PILC REEL
M3696-08	869141B-3B-END-1-2	21	1.00	2.00	2.00	100.0	PILC REEL
M3696-09	869141B-3B-END-1-3	22	1.00	2.00	2.00	100.0	PILC REEL
M3696-10	333000-END-1-1	23	1.00	2.00	2.00	100.0	PILC REEL
M3696-11	333000-END-1-2	24	1.00	2.00	2.00	100.0	PILC REEL
M3696-12	333000-END-2-1	25	1.00	2.00	2.00	100.0	PILC REEL
M3696-13	333000-END-2-2	26	1.00	2.00	2.00	100.0	PILC REEL
M3696-14	333000-END-3-1	27	1.00	2.00	2.00	100.0	PILC REEL
M3696-15	333000-END-3-2	28	1.00	2.00	2.00	100.0	PILC REEL
M3696-16	333000-END-4-1	29	1.00	2.00	2.00	100.0	PILC REEL
M3698-01	FDH095Y-END-1-1	30	1.00	2.00	2.00	100.0	PILC REEL



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 9/10/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 17:20
In Date: 09/09/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:17
Out Date: 09/10/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB116213

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M3698-02	FDH095Y-END-1-2	31	1.00	2.00	2.00	100.0	PILC REEL
M3698-03	183621-END-1-1	32	1.00	2.00	2.00	100.0	PILC REEL
M3698-04	183621-END-1-2	33	1.00	2.00	2.00	100.0	PILC REEL
M3698-05	AFL-05120-END-1-1	34	1.00	2.00	2.00	100.0	PILC REEL
M3698-06	AFL-05120-END-1-2	35	1.00	2.00	2.00	100.0	PILC REEL
M3698-07	AFL-05120-END-2-1	36	1.00	2.00	2.00	100.0	PILC REEL
M3698-08	AFL-05120-END-2-2	37	1.00	2.00	2.00	100.0	PILC REEL

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$