



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 9/16/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 16:30
In Date: 09/15/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:24
Out Date: 09/16/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB116320

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M3733-01	SUP-UST-01-(5-7)	1	1.15	9.96	9.19	91.3	
M3733-02	SUP-UST-01-(10-12)	2	1.15	9.75	8.65	87.2	
M3733-03	SUP-UST-02-(5-7)	3	1.16	9.97	8.75	86.2	
M3733-04	SUP-UST-02-(7-10)	4	1.15	9.93	8.39	82.5	
M3733-05	NAPL-AREA-01-(18-20)	5	1.16	9.72	8.9	90.4	
M3733-06	NAPL-AREA-01-(25-27)	6	1.16	9.97	8.23	80.2	
M3733-07	NAPL-AREA-02-(5-8)	7	1.16	9.96	8.87	87.6	
M3733-08	NAPL-AREA-02-(25-27)	8	1.12	9.83	7.74	76.0	
M3733-09	NAPL-01-(8-10)-END-POINT	9	1.15	9.98	9.24	91.6	
M3733-10	NAPL-02-(8-10)-END-POINT	10	1.15	9.93	9.21	91.8	
M3738-01	TP-16-HEX-CR	11	1.15	9.59	7.4	74.1	
M3738-02	TP-15-HEX-CR	12	1.16	9.97	8.00	77.6	
M3741-01	S-1	13	1.17	9.63	8.88	91.1	
M3741-02	S-2	14	1.12	9.56	8.62	88.9	
M3741-03	S-3	15	1.16	9.96	9.3	92.5	
M3741-04	S-4	16	1.14	9.9	8.41	83.0	
M3741-05	S-5	17	1.16	9.97	7.79	75.3	
M3741-06	S-6	18	1.11	9.97	7.73	74.7	
M3741-07	S-7	19	1.14	9.9	9.27	92.8	
M3741-07D	S-7DUP	20	1.14	9.92	9.29	92.8	
M3743-01	01-A-01-B-01-C	21	1.00	2.00	2.00	100.0	CAULKING SAMPLE
M3743-02	02-A-02-B-02-C	22	1.00	2.00	2.00	100.0	CAULKING SAMPLE
M3743-03	03-A-03-B-03-C	23	1.00	2.00	2.00	100.0	CAULKING SAMPLE
M3743-04	04-A-04-B-04-C	24	1.00	2.00	2.00	100.0	CAULKING SAMPLE
M3743-05	05-A-05-B-05-C	25	1.00	2.00	2.00	100.0	CAULKING SAMPLE
M3743-06	06-A-06-B-06C	26	1.00	2.00	2.00	100.0	CAULKING SAMPLE
M3743-07	07-A-07-B-07-C	27	1.00	2.00	2.00	100.0	CAULKING SAMPLE



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 9/16/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 16:30
In Date: 09/15/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:24
Out Date: 09/16/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB116320

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M3743-08	08-A-08-B-08-C	28	1.00	2.00	2.00	100.0	CAULKING SAMPLE
M3743-09	09-A-09-B-09-C	29	1.00	2.00	2.00	100.0	CAULKING SAMPLE
M3743-10	10-A-10-B-10-C	30	1.00	2.00	2.00	100.0	CAULKING SAMPLE
M3743-11	11-A-11-B-11-C	31	1.00	2.00	2.00	100.0	CAULKING SAMPLE
M3745-01	BC260212LOK-END-1-1	40	1.00	2.00	2.00	100.0	PILC REEL
M3745-02	BC260212LOK-END-1-2	41	1.00	2.00	2.00	100.0	PILC REEL
M3750-01	NAPL-AREA-03- (8-10)	32	1.13	9.68	8.48	86.0	
M3750-02	NAPL-AREA-03- (25-27)	33	1.13	9.59	8.00	81.2	
M3750-03	NAPL-AREA-04- (5-7)	34	1.16	9.72	8.66	87.6	
M3750-04	NAPL-AREA-04- (25-27)	35	1.15	9.95	7.75	75.0	
M3750-05	NAPL-AREA-03- (10-12) -END-POINT	36	1.13	9.6	8.58	88.0	
M3750-06	NAPL-AREA-04- (8-10) -END-POINT	37	1.12	9.78	8.51	85.3	
M3757-01	AU-06-091521	42	1.12	9.69	8.58	87.0	
M3759-01	41095	38	1.00	2.00	2.00	100.0	OIL SAMPLE
M3759-03	20859	39	1.00	2.00	2.00	100.0	OIL SAMPLE

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$