



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 9/21/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 15:15
In Date: 09/20/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:12
Out Date: 09/21/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB116387

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
M3806-01	EWX38	1	1.18	9.68	8.17	82.2	
M3806-02	EWX39	2	1.16	9.89	8.47	83.7	
M3806-03	EWX40	3	1.15	9.93	8.66	85.5	
M3806-04	EWX41	4	1.15	9.94	8.88	87.9	
M3806-05	EWX42	5	1.13	9.52	8.17	83.9	
M3806-06	EWX43	6	1.11	9.67	8.55	86.9	
M3806-07	EWX44	7	1.13	9.82	8.51	84.9	
M3806-08	EWX45	8	1.13	9.69	8.42	85.2	
M3806-09	EWX46	9	1.17	9.81	8.91	89.6	
M3806-10	EWX47	10	1.14	9.63	8.12	82.2	
M3806-11	EWX48	11	1.16	9.81	8.39	83.6	
M3806-12	EWX49	12	1.15	9.51	7.98	81.7	
M3806-13	EWX50	13	1.18	9.78	8.61	86.4	
M3806-14	EWX51	14	1.12	9.64	8.47	86.3	
M3806-15	EWX52	15	1.17	9.9	8.59	85.0	
M3806-16	EWX53	16	1.16	9.9	8.62	85.4	
M3806-17	EWX54	17	1.15	9.95	8.72	86.0	
M3806-18	EWX55	18	1.14	9.94	8.52	83.9	
M3806-19	EWX56	19	1.18	9.83	8.25	81.7	
M3806-20	EWX60	20	1.11	9.87	8.78	87.6	
M3806-21	EWX60D	21	1.11	9.87	8.78	87.6	
M3806-22	EWX60S	22	1.11	9.87	8.78	87.6	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$