



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 9/21/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 15:40
In Date: 09/20/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:15
Out Date: 09/21/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB116389

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M3807-01	EWX57	1	1.14	9.54	8.61	88.9	
M3807-02	EWX58	2	1.12	9.8	8.87	89.3	
M3807-03	EWX59	3	1.14	9.61	8.69	89.1	
M3807-04	EWX61	4	1.16	9.71	8.66	87.7	
M3807-05	EWX62	5	1.19	9.97	8.94	88.3	
M3807-06	EWX63	6	1.13	9.69	8.57	86.9	
M3807-07	EWX64	7	1.16	9.81	9.06	91.3	
M3807-08	EWX65	8	1.17	9.92	9.17	91.4	
M3807-09	EWX66	9	1.16	9.66	8.77	89.5	
M3807-10	EWX67	10	1.16	9.63	8.55	87.2	
M3807-11	EWX68	11	1.15	9.57	8.39	86.0	
M3807-12	EWX69	12	1.13	9.85	8.63	86.0	
M3807-13	EWX70	13	1.14	9.64	8.63	88.1	
M3807-14	EWX71	14	1.14	9.94	8.85	87.6	
M3807-15	EWX72	15	1.16	9.66	8.46	85.9	
M3807-16	EWX73	16	1.11	9.7	8.83	89.9	
M3807-17	EWX74	17	1.14	9.68	8.83	90.0	
M3807-18	EWX75	18	1.1	9.87	8.99	90.0	
M3807-19	EWX76	19	1.18	9.67	8.7	88.6	
M3807-20	EWX84	20	1.16	9.81	8.3	82.5	
M3807-21	EWX84D	21	1.16	9.81	8.3	82.5	
M3807-22	EWX84S	22	1.16	9.81	8.3	82.5	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$