



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 9/23/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 14:00
In Date: 09/22/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:05
Out Date: 09/23/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB116437

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
M3843-01	EWXZ2	1	1.14	9.52	8.38	86.4	
M3843-02	EWXZ3	2	1.17	9.53	8.45	87.1	
M3843-03	EWXZ4	3	1.14	9.54	8.43	86.8	
M3843-04	EWXZ5	4	1.17	9.97	8.97	88.6	
M3843-05	EWXZ6	5	1.13	9.57	8.56	88.0	
M3843-06	EWXZ7	6	1.19	9.7	8.15	81.8	
M3843-07	EWXZ8	7	1.16	9.98	8.53	83.6	
M3843-08	EWXZ9	8	1.18	9.68	8.26	83.3	
M3843-09	EWY00	9	1.13	9.62	8.48	86.6	
M3843-10	EWY01	10	1.19	9.75	8.36	83.8	
M3843-11	EWY02	11	1.14	9.96	8.35	81.7	
M3843-12	EWY03	12	1.17	9.81	8.29	82.4	
M3843-13	EWY03D	13	1.17	9.81	8.29	82.4	
M3843-14	EWY03S	14	1.17	9.81	8.29	82.4	
M3843-15	EWY04	15	1.19	9.69	8.32	83.9	
M3843-16	EWY05	16	1.17	9.7	8.35	84.2	
M3843-17	EWY06	17	1.14	9.62	8.17	82.9	
M3843-18	EWY07	18	1.12	9.5	8.29	85.6	
M3843-19	EWY08	19	1.18	9.75	8.24	82.4	
M3843-20	EWY09	20	1.18	9.97	8.8	86.7	
M3843-21	EWY10	21	1.17	9.55	8.12	82.9	
M3843-22	EWY11	22	1.14	9.61	8.33	84.9	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$