



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 9/23/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 09/22/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:19
Out Date: 09/23/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB116439

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M3861-01	MW-22-20210921-1-2	1	1.14	9.75	8.74	88.3	
M3861-02	MW-22-20210921-2-3	2	1.14	9.82	8.65	86.5	
M3861-03	MW-23-20210921-0-1-WC	3	1.18	9.83	9.22	92.9	
M3861-04	MW-23-20210921-1-2-WC	4	1.17	9.84	8.87	88.8	
M3861-05	MW-23-20210921-2-3-WC	5	1.13	9.56	8.24	84.3	
M3861-06	MW-23-20210921-3-4-WC	6	1.15	9.53	8.34	85.8	
M3861-07	MW-23-20210921-4-5-WC	7	1.13	9.89	8.41	83.1	
M3861-08	MW-23-20210921-0-1	8	1.16	9.79	9.21	93.3	
M3861-09	MW-23-20210921-4-5	9	1.16	9.87	8.02	78.8	
M3862-01	01-A-01-B-01-C	10	1.00	2.00	2.00	100.0	CAULKING SAMPLE
M3862-02	02-A-02-B-02-C	11	1.00	2.00	2.00	100.0	CAULKING SAMPLE
M3862-03	03-A-03-B-03-C	12	1.00	2.00	2.00	100.0	CAULKING SAMPLE
M3862-04	04-A-04-B-04-C	13	1.00	2.00	2.00	100.0	CAULKING SAMPLE
M3862-05	05-A-05-B-05-C	14	1.00	2.00	2.00	100.0	CAULKING SAMPLE
M3862-06	06-A-06-B-06-C	15	1.00	2.00	2.00	100.0	CAULKING SAMPLE
M3867-01	SB-1-(8.0-9.0)	23	1.17	9.78	9.51	96.9	
M3867-02	SB-2-(8.5-9.5)	24	1.19	9.76	9.32	94.9	
M3867-03	SB-3-(7.0-8.0)	25	1.18	9.82	9.35	94.6	
M3867-04	SB-4-(5.0-6.0)	26	1.14	9.96	9.53	95.1	
M3867-05	SB-5-(1.5-1.5)	27	1.13	9.71	9.15	93.5	
M3867-05D	SB-5-(1.5-1.5) DUP	28	1.16	9.73	9.17	93.5	
M3868-01	JC-509	16	1.00	2.00	2.00	100.0	WIPE SAMPLE
M3869-01	SUM-9-22-21	17	1.00	2.00	2.00	100.0	WIPE SAMPLE
M3870-01	TP-C	18	1.14	9.86	8.98	89.9	
M3870-03	TP-C	19	1.14	9.84	9.06	91.0	
M3871-01	SS01-BL-(0.0)	20	1.12	9.62	9.00	92.7	
M3871-02	SS02-BL-(0.0)	21	1.1	9.87	9.04	90.5	
M3871-03	SS03-BL-(0.0)	22	1.1	9.71	8.99	91.6	



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 9/23/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 17:00
In Date: 09/22/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:19
Out Date: 09/23/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB116439

Lab ID	Client SampleID	Dish #	Dish Wt (g) (A)	Dish + Sample Wt (g) (B)	Dish+Dry Sample Wt (g) (C)	% Solid	Comments
M3872-01	SP-LPA-01	29	1.15	9.96	8.03	78.1	
M3872-03	SP-LPA-01	30	1.15	9.96	8.03	78.1	
M3875-01	922A	37	1.00	2.00	2.00	100.0	WIPE SAMPLE
M3875-02	922B	38	1.00	2.00	2.00	100.0	WIPE SAMPLE
M3875-03	20861	39	1.00	2.00	2.00	100.0	OIL SAMPLE
M3875-05	20862	40	1.00	2.00	2.00	100.0	OIL SAMPLE
M3875-07	20863	41	1.00	2.00	2.00	100.0	OIL SAMPLE
M3875-09	20860	42	1.00	2.00	2.00	100.0	OIL SAMPLE
M3875-11	20883	52	1.13	9.98	3.46	26.3	SLUDGE SAMPLE
M3875-13	40441	53	1.00	2.00	2.00	100.0	DEBRIS
M3881-01	OK-01-092221	31	1.18	9.82	8.94	89.8	
M3881-02	OK-01-092221-EPH-2	32	1.17	9.84	8.9	89.2	
M3886-01	COMP-1	33	1.12	9.85	8.72	87.1	
M3886-02	COMP-2	34	1.13	9.67	8.05	81.0	
M3886-02D	COMP-2DUP	35	1.16	9.69	8.07	81.0	
M3887-01	SVN-WRAP-01	36	1.00	2.00	2.00	100.0	wrapping
M3888-01	922C	43	1.00	2.00	2.00	100.0	WIPE SAMPLE
M3888-02	922D	44	1.00	2.00	2.00	100.0	WIPE SAMPLE
M3889-01	SUP-UST-03-(5-7)	45	1.15	9.98	9.19	91.1	
M3889-02	SUP-UST-03-(8-10)	46	1.13	9.93	8.73	86.4	
M3889-03	SUP-UST-04-(6-8)	47	1.13	9.57	8.15	83.2	
M3889-04	SUP-UST-04-(13-15)	48	1.13	9.78	7.08	68.8	
M3891-01	BW-12-20210922-0-1-WC	49	1.16	9.63	8.88	91.1	
M3891-02	BW-12-20210922-1-2-WC	50	1.16	9.66	8.91	91.2	
M3891-03	BW-12-20210922-2-3-WC	51	1.13	9.74	9.02	91.6	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$