



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 9/24/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 13:25
In Date: 09/23/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:45
Out Date: 09/24/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID:

QC:LB116458

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
M3813-01	GB6S7	1	1.12	4.69	4.3	89.1	
M3813-02	GB6S8	2	1.16	9.7	8.52	86.2	
M3813-03	GB6S9	3	1.16	9.93	9.16	91.2	
M3813-04	GB6T0	4	1.16	9.91	8.74	86.6	
M3813-05	GB6T0MS	5	1.16	9.91	8.74	86.6	
M3813-06	GB6T0MSD	6	1.16	9.91	8.74	86.6	
M3813-07	GB6T1	7	1.13	9.93	8.16	79.9	
M3813-08	GB6T2	8	1.16	9.97	8.62	84.7	
M3813-09	GB6T3	9	1.18	9.87	8.3	81.9	
M3813-10	GB6T4	10	1.12	9.54	8.05	82.3	
M3813-11	GB6T5	11	1.16	9.97	8.56	84.0	
M3813-12	GB6T6	12	1.14	9.93	8.91	88.4	
M3813-13	GB6T7	13	1.16	9.98	8.76	86.2	
M3813-14	GB6T8	14	1.15	9.8	8.12	80.6	
M3813-15	GB6T9	15	1.1	9.76	8.59	86.5	
M3813-16	GB6W0	16	1.18	9.97	9.36	93.1	
M3813-17	GB6W1	17	1.13	9.82	8.8	88.3	
M3813-18	GB6W2	18	1.00	2.00	2.00	100.0	P.E. SAMPLE
M3813-19	VHBLK001	19	1.00	2.00	2.00	100.0	VHBLK

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$