



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 9/24/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 13:30
In Date: 09/23/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:47
Out Date: 09/24/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB116459

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
M3837-01	EW7T1	1	1.18	9.6	7.89	79.7	
M3837-02	EW7T2	2	1.15	9.95	8.5	83.5	
M3837-03	EW7T3	3	1.18	9.61	7.89	79.6	
M3837-04	EW7T4	4	1.18	9.74	7.42	72.9	
M3837-05	EW7T5	5	1.19	9.74	7.85	77.9	
M3837-06	EW7T6	6	1.14	9.53	7.88	80.3	
M3837-07	EW7W1	7	1.16	9.84	8.24	81.6	
M3837-08	EW7W2	8	1.16	9.72	8.64	87.4	
M3837-09	EW7W3	9	1.18	9.71	7.97	79.6	
M3837-10	EW7W4	10	1.16	9.87	8.00	78.5	
M3837-11	EW7W5	11	1.15	9.69	8.85	90.2	
M3837-12	EW7W6	12	1.17	9.65	9.04	92.8	
M3837-13	EW7W7	13	1.14	9.98	8.8	86.7	
M3837-14	EW7W8	14	1.18	9.83	8.19	81.0	
M3837-15	EW7W9	15	1.14	9.69	8.07	81.1	
M3837-16	VHBLK001	16	1.00	2.00	2.00	100.0	VHBLK
M3837-17	EW7W1MS	17	1.16	9.84	8.24	81.6	
M3837-18	EW7W1MSD	18	1.16	9.84	8.24	81.6	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$