



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 9/24/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 13:25
In Date: 09/23/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:52
Out Date: 09/24/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB116464

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
M3847-02	EX2W9	1	1.16	9.61	7.91	79.9	
M3847-03	EX2X0	2	1.12	9.83	8.92	89.6	
M3847-04	EX2X1	3	1.14	9.86	8.42	83.5	
M3847-05	EX2X2	4	1.11	9.82	8.56	85.5	
M3847-06	EX2X3	5	1.15	9.97	8.87	87.5	
M3847-07	EX2X4	6	1.12	9.64	8.44	85.9	
M3847-08	EX2X5	7	1.12	9.98	8.64	84.9	
M3847-09	EX2X6	8	1.16	9.71	8.3	83.5	
M3847-10	EX2X7	9	1.19	9.63	8.26	83.8	
M3847-11	EX2X8	10	1.18	9.78	8.47	84.8	
M3847-12	EX2X9	11	1.16	9.57	8.15	83.1	
M3847-13	EX2Y0	12	1.15	9.94	9.46	94.5	
M3847-14	EX2Y1	13	1.13	9.8	9.05	91.3	
M3847-15	EX2Y2	14	1.13	9.52	8.32	85.7	
M3847-16	EX2Y3	15	1.17	9.77	8.77	88.4	
M3847-17	EX2Y4	16	1.1	9.89	8.82	87.8	
M3847-18	EX2Y5	17	1.12	9.87	8.88	88.7	
M3847-19	EX2Y6	18	1.12	9.86	8.89	88.9	
M3847-20	EX2Z0	19	1.17	9.85	8.54	84.9	
M3847-21	EX2Z0D	20	1.17	9.85	8.54	84.9	
M3847-22	EX2Z0S	21	1.17	9.85	8.54	84.9	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$