



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 9/24/2021

OVENTEMP IN Celsius(°C): 107
Time IN: 13:40
In Date: 09/23/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:10
Out Date: 09/24/2021
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-4
Thermometer ID: %SOLIDS-OVEN

QC:LB116466

Lab ID	Client SampleID	Dish #	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish+Dry Sample Wt(g) (C)	% Solid	Comments
M3849-01	EX308	1	1.17	9.83	8.99	90.3	
M3849-02	EX309	2	1.12	9.8	8.93	90.0	
M3849-03	EX310	3	1.18	9.72	8.82	89.5	
M3849-04	EX311	4	1.19	9.63	8.69	88.9	
M3849-05	EX312	5	1.14	9.63	9.06	93.3	
M3849-06	EX313	6	1.14	9.71	9.09	92.8	
M3849-07	EX314	7	1.18	9.67	9.06	92.8	
M3849-09	EX316	8	1.17	9.67	8.85	90.4	
M3849-10	EX317	9	1.19	9.65	8.46	85.9	
M3849-11	EX318	10	1.18	9.91	8.77	86.9	
M3849-12	EX319	11	1.16	9.51	8.33	85.9	
M3849-13	EX320	12	1.15	9.62	8.89	91.4	
M3849-14	EX321	13	1.12	9.83	9.06	91.2	
M3849-15	EX322	14	1.18	9.75	9.00	91.2	
M3849-16	EX323	15	1.14	9.7	8.87	90.3	
M3849-17	EX324	16	1.19	9.66	8.52	86.5	
M3849-18	EX325	17	1.15	9.96	8.86	87.5	
M3849-19	EX326	18	1.17	9.72	8.78	89.0	
M3849-20	EX338	19	1.14	9.94	8.38	82.3	
M3849-21	EX338D	20	1.14	9.94	8.38	82.3	
M3849-22	EX338S	21	1.14	9.94	8.38	82.3	

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$